

100mM dNTP Mix

(Contains 25 mM of each dNTPs)

Catalog No.	Concentration	Size
NT03-10	25 mM of dATP, dCTP, dGTP and dTTP	10x 1 ml
NT03-50	25 mM of dATP, dCTP, dGTP and dTTP	50 ml
NT03-250	25 mM of dATP, dCTP, dGTP and dTTP	250 ml

Description: PCR Nucleotide Mix is a premixed solution containing the sodium salts of dATP, dCTP, dGTP and dTTP, each at a concentration of 10 mM in purified water; the total concentration of nucleotides is 100mM (pH 7.5). This solution is ready to use and is optimized for standard polymerase chain reactions, DNA labeling and sequencing processes.

Applications

- Standard and specialty PCR amplifications
- cDNA synthesis
- qPCR
- Microarrays
- DNA sequencing
- Labeling reactions

Physical Purity: The individual dNTPs contained in the dNTP mix each contain $\geq 99.0\%$ triphosphate as measured by HPLC.

Molecular Weight

dATP: 557.2 g/mol dCTP: 533.1 g/mol
dGTP: 573.1 g/mol dTTP: 548.1 g/mol

Storage Conditions

- Store at -20°C .
- The shelf life is one year from date of receipt under proper storage conditions.
- Avoid freeze/thaw cycles or exposure to frequent temperature changes. It is recommended to aliquot for long term storage
- Mix well prior to use.

Shipping Condition: Blue ice (4°C)

dNTP Mix Reaction Guideline

Reaction Volume	dNTP mix	Final Concentration
50 μl	1 μl	200 μM of each dNTP
25 μl	0.5 μl	200 μM of each dNTP

Standard PCR Composition

Component	Volumes	Final Concentration
10X Reaction Buffer (with 15mM MgCl_2)	5 μl	1X
40 mM Nucleotide Mix (10mM each dNTP)	1 μl	*800 μM
Upstream primer, 20 μM	0.25–2.5 μl	0.1–1 μM
Downstream primer, 20 μM	0.25–2.5 μl	0.1–1 μM
Taq DNA Polymerase, 5U/ μl	0.25-0.50 μl	0.025-0.05 U/ μl
Template DNA	X μl	<250ng
Nuclease-Free Water to final volume of	50 μl	

*The final concentration of each dNTP is 200 μM .



Bioland Scientific LLC, 16612 Maurice Cir., Cerritos, CA 90703

Tel: (562)377-2668 Fax: (562)733-6008 Email: service@bioland-sci.com Online: www.bioland-sci.com