

EasyPrep™ Plasmid Miniprep-II

Manual

Catalog#: PD02-01, PD02-02



For research use only

(March 2018)

Table of Contents

Introduction.....	3
Important Notes.....	3
Before Starting.....	4
Kit Contents	5
Storage and Stability.....	5
Plasmid Miniprep-II Centrifuge Method.....	6
Plasmid Miniprep-II Vacuum Method.....	8
Trouble Shooting Guide	10

Introduction

The EasyPrep™ Plasmid Miniprep-II Kit provides a simple, fast and cost effective method to purify plasmid DNAs from bacterial culture, especially useful for low-copy-number plasmid miniprep or applications that amount of DNA needed is about 30-50 µg. The kit applies DNA-binding-specific matrix, which allows the high efficient binding of DNA, while proteins and other cellular components are easily removed during washing. Plasmid DNAs are easily eluted with Tris buffer or pure water.

The yield of the high-copy plasmid from this kit is about 30-50µg. Plasmid DNAs purified from this kit can be used directly for downstream applications such as transfection, cloning, sequencing, as well as *in vivo* applications.

Important Notes

The yield of plasmid DNA using the Miniprep-II kit depends on:

- 1) The origin of the replication. Samples of high-copy-number and low-copy-number plasmids are shown in the table. The protocols are optimized for high-copy-number plasmid purification. For low-copy-number plasmids, both the culture volume and the buffer volume need to be scaled up 3 to 5 times or use our EasyPrep™ plasmid Miniprep-II kits.
- 2) The size of the plasmids. In general, larger plasmids have higher yields than smaller plasmids if the copy-numbers are similar. For example, the yield of pcDNA3.1 (5.4kb) > pBluescript (3kb) > pUC18/19 (2.7kb).

Plasmid	Origin	High copy	Low copy
pACYC	P15A		10-12
pSC101	pSC101		5
pSuperCos	pMA2		10-20
pBR322	pMA2		15-20
pUC	Muted pMA2	500-700	
pGEM ^R	Muted pMA2	300-400	
pBluescript ^R	ColE1	300-500	

- 3) Culture condition. For maximal yields, the culture should be in late log or early stationary phase of bacteria growth. If the culture overgrows, you will be harvesting more dead bacteria than live cells, which leads to genomic DNA contamination. If the culture is undergrown, then of course, yields are lower than expected. In practice:

- Start culture with fresh colonies, no more than a few days old.
- Add antibiotics to the LB media just before the culture.
- Shaker speed, 250 rpm for 5-12 ml LB medium in 50 ml culture tube or flasks.
- Culture at 37°C for 12-16 hours.
- OD₆₀₀ of the culture should be 2.5-3.5.

Before Starting

Prepare all components and get all necessary materials ready by examining this manual and become familiar with each step.

- **RNase A:** Spin down RNase A vial briefly. Add the whole tube of RNase A to Buffer A1 and mix well before use. Store Buffer A1 at 4°C after use. RNase A final concentration ≥ 100 µg/ml.
- **Buffer A2** could precipitate below room temperature. It is critical to warm up the buffer at 37°C to dissolve the precipitates before use.
- Keep Buffer A2 bottle cap tightly closed to avoid acidification.
- **Add ethanol to DNA Washing Buffer bottle before use:**
Cat#: PD02-01: Add 80 mL 96-100% ethanol and mix
Cat#: PD02-02: Add 216 mL 96-100% ethanol and mix
- Carry out all centrifugations at **room temperature**.

Materials supplied by users:

- 96-100% ethanol
- 1.5 ml microcentrifuge tubes
- High speed microcentrifuge

Kit Contents

Catalog#	PD02-01	PD02-02
Miniprep-II DNA Columns	50	250
Buffer A1	30 ml	130 ml
Buffer A2	30 ml	130 ml
Buffer N1	45 ml	180 ml
DNA Wash Buffer Concentrate*	20 ml	54 ml
Elution Buffer**	10 ml	30 ml
Column Buffer	10 ml	30 ml
RNase A	110 µl	440 µl
2ml Collection tubes	50	250
Manual	1	1

* Add ethanol to DNA Washing Buffer bottle (see Page 4) before use

**Elution Buffer: 10mM Tris.HCl, 1mM EDTA, pH8.5 (TE buffer, pH8.5)

Storage and Stability

Buffer A1 should be stored at 4°C once RNase A is added. All other components can be stored at room temperature. Shelf life: 12 months from the date of purchase.

EasyPrep™ Plasmid Miniprep-II Protocol

(Centrifuge Method)

1. **Harvest 5-12 ml overnight bacterial culture by centrifugation at 5,000 g for 10 min (or 3,500 rpm for 20 min) in a tabletop centrifuge. Pour off the supernatant and blot the inverted tube on a paper towel to remove excess media.**
2. **Add 450 µl Buffer A1 and completely resuspend bacterial pellet by vortexing or pipetting.** *Ensure that RNase A is added into Buffer A1 before use. Complete resuspension is critical for optimal yields.*
3. **Add 450 µl Buffer A2, mix gently but thoroughly by inverting 5-6 times, incubate for approximately 2-4 min until the solution become slightly clear.** *Do not incubate longer than 5 min. Over-incubation causes genomic DNA contamination and plasmid damage. Buffer A2 may form precipitates below room temperature. Warm up at 37°C water-bath to dissolve the precipitates before use.*
4. **Add 600 µl Buffer N1, mix immediately and thoroughly by inverting 5-6 times. Don't vortex!** *It is critical to mix the solution well, if it appears conglobated, more mix is required.*
5. **Centrifuge at 13,000 rpm for 5 min at room temperature.** *If the lysate doesn't appear clear, reverse the tube angle, centrifuge for 5 min and then transfer the lysate to the DNA column.*
6. **Column Equilibration: during centrifuge, add 100 µl Column Buffer to the column, incubate for 2 min, centrifuge at 11,000 rpm for 30-60s.** *This step ensures optimal binding capacity.*
7. **Carefully transfer up to 750µl clear lysate into a DNA column with the collection tube** (do not bring any floating material into the tube). **Centrifuge at 11,000 rpm for 30-60s. Discard the flow-through.** *When using end+ bacteria strains, such as JM series and HA201 derivatives, it's important removes trace nuclease. Add 400 µl PB Buffer into the DNA column, centrifuge at 11,000 rpm for 30-60 s. Discard the flow-through. End- strains such as DH5α, Top10 and XL-1 blue do not require this treatment.*

7. Add **750 µl DNA Wash Buffer** into the DNA column, centrifuge at **11,000 rpm for 30-60s**. Discard the flow-through. *Ensure that the 96-100% ethanol has been added to DNA Wash Buffer according to instructions before use.*
8. Centrifuge at **11,000 rpm for additional 1 min**. *It is critical to remove residual ethanol for optimal elution. Residual ethanol may interfere subsequent applications.*
9. Place the DNA column in a clean **1.5 ml microcentrifuge tube** (Be careful not to transfer any of the ethanol with the DNA column). Add **50-100 µl Elution Buffer** (TE buffer, pH8.5) to the center of the column, incubate for **1 min**, centrifuge at **11,000 rpm for 1 min** to elute the plasmid DNA. *Elution buffer can be replaced with 10mM Tris.HCl, pH8.5 if EDTA interferes subsequent enzyme digestions. Proper pH is essential.*

DNA storage:

Plasmid DNA is stable in TE Buffer if stored at -20°C or below.

Yield and quality of DNA

Determine the DNA absorbance of an appropriate dilution (20 to 50 fold) of the DNA at 260 nm and 280 nm. The DNA concentration is calculated as below:

DNA concentration = $A_{260} \times 50 \times \text{dilution factor}$ (µg/ml)

EasyPrep™ Plasmid Miniprep-II Protocol

(Vacuum Method)

1. Harvest **5-12 ml** overnight bacterial culture by centrifugation at **5,000 g for 10 min** (or **3,500 rpm for 20 min**) in a tabletop centrifuge. Pour off the supernatant and blot the inverted tube on a paper towel to remove excess media.
2. Add **450 µl Buffer A1** and completely resuspend bacterial pellet by vortexing or pipetting. *Ensure that RNase A is added into Buffer A1 before use. Complete resuspension is critical for optimal yields.*
3. Add **450 µl Buffer A2**, mix gently but thoroughly by inverting **5-6 times**, incubate for approximately **2-4 min** until the solution become slightly clear. *Do not incubate longer than 5 min. Over-incubation causes genomic DNA contamination and plasmid damage. Buffer A2 may form precipitates below room temperature. Warm up at 37°C water-bath to dissolve the precipitates before use.*
4. Add **600 µl Buffer N1**, mix immediately and thoroughly by inverting **5-6 times**. Don't vortex! *It is critical to mix the solution well, if it appears conglobated, more mix is required.*
5. Centrifuge at **13,000 rpm for 5 min at room temperature**. *If the lysate doesn't appear clear, reverse the tube angle, centrifuge for 5 min and then transfer the lysate to the DNA column.*
6. During centrifugation, prepare the vacuum manifold as instructed by manufacturer. Insert Miniprep-II columns into the luer connector of the vacuum manifold. Add **100 µl Column Buffer** to the column, incubate for **2 min**, switch on vacuum to draw the solution through the columns.
7. Carefully transfer up to **750 µl clear lysates into Miniprep-II DNA columns** (Do not bring any floating material into the tube). Switch on vacuum to draw the solution through the columns completely, and then switch off vacuum.

Note: when using end+ bacteria strains, such as JM series and HA201 derivatives, it's important to remove trace nuclease. Add 400 µl PB Buffer into the DNA column, centrifuge at 11,000 rpm for 30-60 s. Discard the flow-through. End⁻ strains such as DH5α,

Top10 and XL-1 blue do not require this step.

8. **Add 750 µl DNA Wash Buffer** into the DNA columns. **Switch on vacuum source to draw the solution through the columns completely, and then switch off vacuum source.** *Ensure that the 96-100% ethanol has been added to DNA Wash Buffer according to instructions before use.*
9. **Transfer the DNA columns to 2ml collection tubes. Centrifuge at 11,000 rpm for 1 min .** *It is critical to remove residual ethanol for optimal elution. Residual ethanol may interfere subsequent applications.*
7. **Place each DNA column in a clean 1.5 ml microcentrifuge tube** (Be careful not to transfer any of the ethanol with the DNA column). **Add 50-100 µl Elution Buffer** (TE buffer, pH8.5) **to the center of the column, incubate for 1 min, centrifuge at 11,000 rpm for 1 min to elute the plasmid DNA.** *Elution buffer can be replaced with 10mM Tris.HCl, pH8.5 if EDTA interferes subsequent enzyme digestions.*

DNA storage:

Plasmid DNA is stable in TE Buffer if stored at -20°C or below.

Yield and quality of DNA

Determine the DNA absorbance of an appropriate dilution (20 to 50 fold) of the DNA at 260 nm and 280 nm. The DNA concentration is calculated as below:

$$\text{DNA concentration} = A_{260} \times 50 \times \text{dilution factor} (\mu\text{g/ml})$$

Trouble Shooting Guide

Problem	Possible Reason	Suggested Improvement
Low Yield	Poor Cell lysis	Resuspend pellet thoroughly by vortexing and pipetting before adding Buffer A2. Make fresh Buffer A2 if the cap had not been closed tightly.
	Bacterial culture over-grown or not fresh	Grow bacterial 12-16 h. If the culture is not used the same da, spin down cultures and store the pellet at -20°C. Do not store culture at 4°C overnight.
	Low copy-number plasmid	Increase culture volume (up to 15 ml for Miniprep-III). Increase the volume of Buffer A1, A2, N1 as instructed in the protocol.
No DNA	Plasmid lost in Host <i>E. coli</i>	Prepare fresh culture from fresh colonies.
Genomic DNA contamination	Over incubation after adding Buffer A2	Do not vortex or mix aggressively after adding Buffer A2. Do not incubate more than 5 min after adding solution A2.
RNA contamination	RNase A not added to Buffer A1 or Buffer A1 not stored at 4°C after adding RNase A.	Add RNase A to Buffer A1: 100 µl RNase A (30mg/ml) for 30 ml Buffer A1.
Plasmid DNA floats out of wells while running in agarose gel, DNA doesn't freeze or smell of ethanol	Residual ethanol not completely removed from column	Centrifuge or vacuum again after washing step. Make sure that no residual ethanol remaining on the membrane before eluting the plasmid DNA.

Limited use and warranty

This product is warranted to perform as described in its labeling and in Bioland's literature when used in accordance with instructions. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Bioland. Bioland's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of Bioland, to replace the products. Bioland shall have no liability for any direct, indirect, consequential, or incidental damage arising out of the use, the results of use, or the inability to use it product.

Bioland Scientific LLC

**14925 Paramount Blvd., Suite C
Paramount, CA 90723
USA**

Tel: (877) 603-8882; (562) 602-8882

Fax: (562) 733-6008

Email: service@bioland-sci.com
order@bioland-sci.com

Visit our web at www.bioland-sci.com and learn more about
Bioland products