

## 2x Multiplex qPCR Mastermix for TaqMan Probe (No ROX)

Cat. No.	Product Description	Size
QPTM01-00	2x Multiplex qPCR Mastermix for TaqMan probe	1 ml
QPTM01-01	2x Multiplex qPCR Mastermix for TaqMan probe	5x1 ml
QPTM01-02	2x Multiplex qPCR Mastermix for TaqMan probe	50 ml

### Intended Use

- The 2x Multiplex qPCR Mastermix is used for real-time qualitative and quantitative multiplex qPCR with TaqMan probes for up to four targets.
- The Mastermix is a premixed, 2X concentrated solution that has all the components except for gene-specific primers, probes and templates

### Kit Characterizations

- The kit is designed for multiplex qPCR with TaqMan probes.
- The kit uses Super *Taq*-Probe DNA polymerase specially engineered for TaqMan probes, which increases 5' to 3' exonuclease efficiency and produces S-shaped curve.
- Two to four pairs of gene-specific primers can be applied in one reaction.
- The concentrations of the primers and probes are variable depending on specific assays and thermocycling protocols (Table 1).
- The preferable PCR product size is ≤150bp.
- The kit has three formulations of No ROX, Low ROX or High ROX concentrations for your choice (see Table 2).

### Transportation and storage

- The kit can be transported at below 4°C for up to 3 days.
- The kit should be kept stable in the dark at -20°C for ≤24 months with ≤10 times of freeze-thaw cycles. The kit can be stored at 4°C for a week.

### Setup Reaction and Thermocycling

- Thaw 2x Mastermix and other reaction components at room temperature, mix each component, centrifuge and then place on ice.
- Set up reactions (Table 1).
- Seal tubes with flat, optically transparent caps or seal plate with optically transparent film.
- Mix and then briefly centrifuge the tubes or plate.
- Program PCR instrument with indicated thermo-cycling protocol.
- Load PCR tubes or plate and start to run.
- Perform data analysis according to the PCR instrument instructions.

**Table 1. Setting up a 20µL or 10µL reaction**

Component	Volume per 20µL	Volume per 10µL	Final concentration
2x Mastermix	10µL	5µL	1X
Primers <sup>a</sup>	Variable	Variable	Each 150-400nM
TaqMan probes <sup>b</sup>	Variable	Variable	Each 150-250nM
DNA templates <sup>c</sup>	Variable	Variable	≤500ng human genomic DNA/20µL
H <sub>2</sub> O	To 20µL	To 10µL	

### Footnotes of Table 1

- Each primer's  $T_m$  should be designed ≥60°C, preferably between 62°C to 65°C, using primer3 software for high efficiency and specificity.
- Each probe's  $T_m$  should be 8-10°C higher than the primer's  $T_m$ , preferably between 70-75°C.
- DNA templates should be extracted by a qualified silica-based kit and eluted with low EDTA TE buffer (10mM Tris-HCl, 0.1mM EDTA, pH 8.0-8.3).



## 2x Multiplex qPCR Mastermix for TaqMan Probe (No ROX)

**Table 2. Compatible instruments**

qPCR Instrument	ROX required by instrument	Passive dye setup
Bio-Rad® iQ™5, CFX96, CFX384, Opticon Roche Lightcycler® Qiagen Rotor-Gene™ Eppendorf Mastercycler® Cepheid® SmartCycler®	Not recommended	Not necessary
Applied Biosystems® 7500, 7500 Fast, QuantStudio™, ViiA7™, Agilent Mx™	Low ROX (50nM final concentration)	Turn on ROX passive reference dye
Applied Biosystems® 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast, StepOne™, StepOnePlus™	High ROX (500nM final concentration)	Turn on ROX passive reference dye

**Table 3. Standard Thermocycling Protocol**

Stage	Temperature	Period	Number of cycles
I	95°C	2min	1
II	95°C	10sec	35-40
	60°C, signal acquisition	60sec	

**Table 4. Three-Step Thermocycling Protocol**

Stage	Temperature	Period	Number of cycles
I	95°C	2min	1
II	95°C	10sec	35-40
	60°C	30sec	
	68-72°C, signal acquisition	30sec	

**Footnotes of Tables 3 and 4**

- The three-step thermocycling protocol in Table 4 increases overall polymerase activity by 50%, a more effective protocol than Table 3.
- The primer concentration used in Tables 3 and 4 is typically 0.15-0.2µM.

### Related Products

- 2x qPCR Mastermix for TaqMan Probe
- 2x Fast qPCR Mastermix-SYBR Green

### Precautions

If you order a “**No ROX**” Mastermix but you have an Applied Biosystems or Thermo Fisher instrument, please **turn off ROX passive reference dye button** when setup assays.