

# Safety Data Sheet

Trade Name: BL Buffer

Revision Date: 4/28/2016

## SECTION 1. IDENTIFICATION OF PRODUCT AND COMPANY:

**Product Name:** BL Buffer

### Manufacturer or Supplier Details:

Company Name: Omega Bio-tek, Inc.  
Address: 400 Pinnacle Way, Suite 450  
Norcross, GA 30071  
USA  
  
Telephone: 1-800-832-8896  
Email: info@omegabiotek.com  
Emergency Telephone: CHEMTREC  
USA & Canada: 1-800-424-9300  
Outside USA & Canada: 1-703-527-3887

### Recommended Use of the Chemical and Restrictions on Use:

Recommended Use: For research use only.

## SECTION 2. HAZARDS IDENTIFICATION

### Hazard Classification:

|                 |   |   |
|-----------------|---|---|
| Health Hazards: | Acute toxicity (oral)<br>Skin corrosion/irritation<br>Serious eye damage/eye irritation<br>Specific target organ toxicity - single exposure | Category 4<br>Category 2<br>Category 2A<br>Category 3 |
|-----------------|---|---|

### GHS Label Element:

Hazard Symbol:



Signal Word:

Warning

### Hazard Statements:

H302 - Harmful if swallowed.  
H315 - Causes skin irritation.  
H319 - Causes serious skin irritation.  
H332 - Harmful if inhaled.

# Safety Data Sheet

Trade Name: BL Buffer

Revision Date: 4/28/2016

## Precautionary Statements:

P261 - Avoid breathing dust/fumes/gas/mist/vapors/spray.  
P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P311 - Call a poison control center or doctor/physician.

## Potential Health Effects:

Aggravated Medical Conditions: May cause nervous system disturbances.  
Symptoms of Overexposure:  
Inhalation: Toxic if inhaled. May cause respiratory tract irritation.  
Skin: Harmful if absorbed through skin. May cause skin irritation.  
Eyes: Causes eye irritation.  
Ingestion: Harmful if swallowed.

## Carcinogenicity:

IARC No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.  
ACGIH No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH.  
OSHA No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.  
NTP No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.

## SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS:

Substance or Mixture: Mixture of the substances listed below with non-hazardous additions.

Synonyms/Common Names: Guanidinium chloride, Guanidinium monohydrochloride, Guanidinium hydrochloride

| Chemical Name                                      | CAS No. | Concentration (%) |
|--|---------|-------------------|
| Guanidine Hydrochloride<br>$\text{CH}_6\text{N}_3$ | 50-01-1 | 50-100%           |

# Safety Data Sheet

Trade Name: BL Buffer

Revision Date: 4/28/2016

## SECTION 4. FIRST AID MEASURES:

|   |  |
|---|--|
| <b>General Advice:</b>  | Exit to a safe area. Consult a physician. Seek medical attention.                                      |
| <b>If Inhaled:</b>  | Move to fresh air. If not breathing, give artificial respiration. Seek medical attention.              |
| <b>In Case of Skin Contact:</b>   | Exit to safe area. Wash thoroughly with soap and water. Seek medical attention.                        |
| <b>In Case of Eye Contact:</b>  | Flush eyes with clean water for a minimum of 15 minutes, keeping eyes open. Seek medical attention.    |
| <b>If Swallowed:</b>  | Rinse mouth with water. Never give anything by mouth to an unconscious person. Seek medical attention. |
| <b>Most Important Symptoms and Effects, Acute and Delayed:</b>                | Irritating to eyes, respiratory system and skin. Narcotic effect.                                      |
| <b>Recommendations for Immediate Medical Attention and Special Treatment:</b> | Treat symptomatically. Symptoms may be delayed.  |

## SECTION 5. FIREFIGHTING MEASURES:

|   |  |
|---|--|
| <b>Conditions of Flammability:</b>                    | Not flammable or combustible.  |
| <b>Suitable Extinguishing Media:</b>                  | Dry powder, CO <sub>2</sub> , water spray or regular foam.                         |
| <b>Special Protective Equipment for Firefig ters:</b> | Wear self-contained breathing apparatus, if necessary.                             |
| <b>Hazardous Combustion Products:</b>                 | Carbon oxides (CO, CO <sub>2</sub> ), nitrogen oxides (NOx), hydrogen chloride gas |
| <b>Further Information:</b>                           | Prevent firefighting water from entering surface water or groundwater.             |

# Safety Data Sheet

Trade Name: BL Buffer

Revision Date: 4/28/2016

## SECTION 6. ACCIDENTAL RELEASE MEASURES:

|   |   |
|---|---|
| <b>Personal Precautions:</b>                                  | Wear respiratory protection. Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. See Section 8: Exposure Controls/Personal Protection for more details. |
| <b>Environmental Precautions:</b>                             | Prevent further leakage if area is safe. Do not let product enter drains. Avoid discharge into the environment.   |
| <b>Emergency Procedures:</b>                                  | Evacuate personnel to safe areas and follow emergency response protocols.   |
| <b>Methods and Materials for Containment and Cleaning Up:</b> | Use inert absorbent material to soak up spill. Dispose of hazardous waste. Do not flush with water.   |

## SECTION 7. HANDLING AND STORAGE:

|                                       |  |
|---------------------------------------|--|
| <b>Precautions for Safe Handling:</b> | Handle in accordance with good chemical hygiene and safety practices. Protect from moisture. Avoid contact with eyes, skin and clothing. Avoid inhalation of vapors or mist. |
| <b>Conditions for Safe Storage:</b>   | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.  |
| <b>Incompatible Products:</b>         | Strong oxidizing agents  |

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

**Appropriate Engineering Controls:** Adequate ventilation and access to eye washing stations are required.

### Recommendations for Personal Protective Measures:

|                                  |   |
|----------------------------------|---|
| <b>Hand Protection:</b>          | Use gloves when handling this product. Inspect gloves prior to use for any visible damage. Use proper glove removal technique to prevent contact with skin. Dispose of gloves in accordance with applicable laws. Wash and dry hands after use. |
| <b>Eye Protection:</b>           | Use a face shield and/or safety goggles when handling this product. Use eyewear that has been approved by NIOSH.  |
| <b>Skin and Body Protection:</b> | Wear appropriate laboratory attire when handling this product, such as gloves, long pants, closed-toe shoes and a laboratory coat.  |

# Safety Data Sheet

Trade Name: BL Buffer

Revision Date: 4/28/2016

|                         |   |
|-------------------------|---|
| Respiratory Protection: | If engineering controls do not maintain airborne concentrations below recommended exposure limits (where applicable) or to an acceptable level (in countries where exposure limits have not been established), an approved respirator must be used. |
| Hygiene Measures:       | Avoid contact with eyes, skin and clothing. Wash hands thoroughly after using this product. Do not eat, drink or smoke when handling this product.  |

**Special Requirements for Personal Protective Equipment, Protective Clothing, Respirators, etc:**

**Components with Workplace Control Parameters:**

| Components | CAS No. | Value | Control Parameters | Basis |
|------------|---------|-------|--------------------|-------|
|------------|---------|-------|--------------------|-------|

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES:

|   |                   |
|---|-------------------|
| <b>Appearance:</b>  |                   |
| Form:   | Liquid            |
| Color:  | Yellow            |
| Odor:   | No data available |
| Odor Threshold:   | No data available |
| pH:   | 3.0-5.0           |
| Melting Point/Freezing Point:                                   | No data available |
| Initial Boiling Point/Boiling Range:                            | No data available |
| Flash Point:  | No data available |
| Evaporation Rate:   | No data available |
| Flammability (Solid, Gas):                                      | No data available |
| <b>Upper/Lower Limits on Flammability and Explosive Limits:</b> |                   |
| Flammability Limit (Upper):                                     | No data available |
| Flammability Limit (Lower):                                     | No data available |
| Explosive Limit (Upper):  | No data available |
| Explosive Limit (Lower):  | No data available |
| Vapor Pressure:   | No data available |
| Vapor Density:  | No data available |

# Safety Data Sheet

Trade Name: BL Buffer

Revision Date: 4/28/2016

|   |                   |
|---|-------------------|
| <b>Relative Density:</b>                        | No data available |
| <b>Solubility:</b>                              |                   |
| Solubility in Water:                            | No data available |
| Solubility (Other):                             | No data available |
| <b>Partition Coefficient (n-octanol/water):</b> | No data available |
| <b>Auto-Ignition Temperature:</b>               | No data available |
| <b>Decomposition Temperature:</b>               | No data available |
| <b>Viscosity:</b>                               | No data available |

## SECTION 10. STABILITY AND REACTIVITY:

|  |   |
|--|---|
| <b>Reactivity:</b>                         | No dangerous reaction known under conditions of normal use.                   |
| <b>Chemical Stability:</b>                 | Stable under recommended storage conditions.                                  |
| <b>Possibility of Hazardous Reactions:</b> | Hazardous polymerization does not occur.                                      |
| <b>Conditions to Avoid:</b>                | Contact with incompatible materials, excessive heat.                          |
| <b>Materials to Avoid:</b>                 | Strong oxidizing agents, strong acids, inorganic salts, chlorinated compounds |
| <b>Hazardous Decomposition Products:</b>   | Carbon oxides, nitrogen oxides (NOx), hydrogen chloride gas                   |

## SECTION 11. TOXICOLOGICAL INFORMATION:

### Information on Likely Routes of Exposure:

|             |   |
|-------------|---|
| Inhalation: | Harmful if inhaled. May cause respiratory tract irritation. |
| Ingestion:  | Harmful if swallowed.                                       |
| Skin:       | Causes skin irritation.                                     |
| Eyes:       | Causes serious eye irritation.                              |

### Information on Toxicological Effects:

|                    |                      |
|--------------------|----------------------|
| Acute Toxicity:    |                      |
| Oral LD50:         | 1,900 mg/kg (Rat)    |
| Inhalation LC50:   | No data available    |
| Dermal LD50:       | 2,000 mg/kg (Rabbit) |
| Other Information: | No data available    |

# Safety Data Sheet

Trade Name: BL Buffer

Revision Date: 4/28/2016

|                                    |  |
|------------------------------------|--|
| Skin Corrosion/Irritation:         | Causes skin irritation.  |
| Serious Eye Damage/Irritation:     | Causes serious eye irritation.   |
| Respiratory or Skin Sensitization: | Components of this product may be absorbed into the body through the skin. Irritating to skin and mucus membranes. |
| Germ Cell Mutagenicity:            | No mutagenic components identified.  |

## **Carcinogenicity:**

|        |  |
|--------|--|
| IARC:  | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.  |
| ACGIH: | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH. |
| NTP:   | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.   |
| OSHA:  | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.  |

**Reproductive Toxicity:** No component toxic to reproduction.

**Teratogenicity:** No data available.

**Specific Target Organ Toxicity - Single Exposure:** Respiratory tract irritation.

**Specific Target Organ Toxicity - Repeated Exposures:** No data available.

**Aspiration Hazard:** No data available.

**Signs and Symptoms of Exposure:** Irritating to eyes, respiratory system and skin. Narcotic effect.

**Synergistic Effects:** No data available.

**Additional Information:** None known.

## **SECTION 12. ECOLOGICAL INFORMATION:**

|                                       |  |
|---------------------------------------|--|
| <b>Toxicity:</b>                      | No data available.   |
| <b>Persistence and Degradability:</b> | No data available.   |
| <b>Bioaccumulative Potential:</b>     | No data available.   |
| <b>Mobility in Soil:</b>              | Product is water soluble and may spread in water systems.  |
| <b>PBT and vPvB Assessment:</b>       | No data available.   |
| <b>Other Adverse Effects:</b>         | The product components are not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment. |

# Safety Data Sheet

Trade Name: BL Buffer

Revision Date: 4/28/2016

## SECTION 13. DISPOSAL CONSIDERATIONS:

**Product:**

Disposal should be in accordance with applicable regional, national and local laws and regulations.

S29 - Do not empty into drains.

S57 - Use appropriate container to avoid environmental contamination.

Refer to Section 7: Handling and Storage and Section 8: Exposure Control/Personal Protection for additional handling information and protection of employees.

**Contaminated Packaging:**

Disposal must be made according to official regulations.

To minimize exposure, please refer to Section 8: Exposure Control/Personal Protection of this SDS.

## SECTION 14. TRANSPORT INFORMATION:

**DOT:**

Not regulated.

**IMDG:**

Not regulated.

**IATA:**

Not regulated.

## SECTION 15. REGULATORY INFORMATION:

**US Federal Regulations:**

OSHA Hazards:

Aggravated Medical  
Conditions:

May cause nervous system disturbances.

Symptoms of  
Overexposure:

Inhalation:

Toxic if inhaled. May cause respiratory tract irritation.

Skin:

Harmful if absorbed through skin. May cause skin irritation.

Eyes:

Causes eye irritation.

Ingestion:

Harmful if swallowed.

SARA 302 Components:

None present or none present in regulated quantities.

SARA 313 Components:

None present or none present in regulated quantities.

SARA 311/312 Hazards:

Acute health hazard, chronic health hazard



# Safety Data Sheet

Trade Name: BL Buffer

Revision Date: 4/28/2016

## US State Regulations:

|                             |  |
|-----------------------------|--|
| California Prop. 65:        | No ingredient regulated by CA Prop 65 present.           |
| New Jersey & Community RTK: | No ingredient regulated by NJ Right to Know Act present. |
| Massachusetts RTK:          | No ingredient regulated by MA Right to Know Act present. |
| Pennsylvania RTK:           | No ingredient regulated by PA Right to Know Act present. |
| Rhode Island RTK:           | No ingredient regulated by RI Right to Know Act present. |

## SECTION 16. OTHER INFORMATION:

|                      |  |
|----------------------|--|
| Issue Date:          | 8/19/2015  |
| Revision Date:       | 4/28/2016  |
| Version:             | 5.0  |
| Further Information: | No data available.   |
| Disclaimer:          | The above information is believed to be correct but does not claim to be all inclusive. This document shall be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Omega Bio-tek, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. |

# Safety Data Sheet

Trade Name: DNA Wash Buffer

Revision Date: 4/29/2016

## SECTION 1. IDENTIFICATION OF PRODUCT AND COMPANY:

**Product Name:** DNA Wash Buffer

### Manufacturer or Supplier Details:

|                      |  |
|----------------------|--|
| Company Name:        | Omega Bio-tek, Inc.  |
| Address:             | 400 Pinnacle Way, Suite 450<br>Norcross, GA 30071<br>USA                         |
| Telephone:           | 1-800-832-8896   |
| Email:               | info@omegabiotek.com   |
| Emergency Telephone: | CHEMTREC<br>USA & Canada: 1-800-424-9300<br>Outside USA & Canada: 1-703-527-3887 |

### Recommended Use of the Chemical and Restrictions on Use:

|                  |                        |
|------------------|------------------------|
| Recommended Use: | For research use only. |
|------------------|------------------------|

## SECTION 2. HAZARDS IDENTIFICATION

### Hazard Classification:

|                 |   |
|-----------------|---|
| Health Hazards: | No known OSHA hazards.<br>Not a dangerous substance according to GHS. |
|-----------------|---|

### GHS Label Element:

|                |     |
|----------------|-----|
| Hazard Symbol: | N/A |
|----------------|-----|

|                    |                                       |
|--------------------|---------------------------------------|
| Hazard Statements: | Not a hazardous substance or mixture. |
|--------------------|---------------------------------------|

## SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS:

|                       |                                      |
|-----------------------|--------------------------------------|
| Substance or Mixture: | Mixture of non-hazardous substances. |
|-----------------------|--------------------------------------|

|                        |     |
|------------------------|-----|
| Synonyms/Common Names: | N/A |
|------------------------|-----|

# Safety Data Sheet

Trade Name: DNA Wash Buffer

Revision Date: 4/29/2016

## SECTION 4. FIRST AID MEASURES:

|   |  |
|---|--|
| <b>General Advice:</b>  | Wear protective gloves/protective clothing/eye protection/face protection. |
| <b>If Inhaled:</b>  | If not breathing, give artificial respiration.                             |
| <b>In Case of Skin Contact:</b>   | Wash with soap and water.  |
| <b>In Case of Eye Contact:</b>  | Remove contact lenses, if present and easy to do. Continue rinsing.        |
| <b>If Swallowed:</b>  | N/A  |
| <b>Most Important Symptoms and Effects, Acute and Delayed:</b>                | N/A  |
| <b>Recommendations for Immediate Medical Attention and Special Treatment:</b> | Treat symptomatically. Symptoms may be delayed.                            |

## SECTION 5. FIREFIGHTING MEASURES:

|   |  |
|---|--|
| <b>Conditions of Flammability:</b>                    | Not flammable or combustible.  |
| <b>Suitable Extinguishing Media:</b>                  | Dry powder, CO <sub>2</sub> , water spray or regular foam.             |
| <b>Special Protective Equipment for Firefighters:</b> | Wear self-contained breathing apparatus, if necessary.                 |
| <b>Hazardous Combustion Products:</b>                 | N/A  |
| <b>Further Information:</b>                           | Prevent firefighting water from entering surface water or groundwater. |

# Safety Data Sheet

Trade Name: DNA Wash Buffer

Revision Date: 4/29/2016

## SECTION 6. ACCIDENTAL RELEASE MEASURES:

|   |   |
|---|---|
| <b>Personal Precautions:</b>                                  | Avoid inhalation of vapors or mist. See Section 8: Exposure Controls/ Personal Protection for more information. |
| <b>Environmental Precautions:</b>                             | N/A   |
| <b>Emergency Procedures:</b>                                  | Evacuate personnel to safe areas and follow emergency response protocols.                                       |
| <b>Methods and Materials for Containment and Cleaning Up:</b> | Use inert absorbent material to soak up spill.  |

## SECTION 7. HANDLING AND STORAGE:

|                                       |  |
|---------------------------------------|--|
| <b>Precautions for Safe Handling:</b> | Handle in accordance with good chemical hygiene and safety practices. Protect from moisture. Avoid contact with eyes, skin and clothing. Avoid inhalation of vapors or mist. |
| <b>Conditions for Safe Storage:</b>   | Keep container tightly closed in a dry, well-ventilated place.   |
| <b>Incompatible Products:</b>         | N/A  |

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

**Appropriate Engineering Controls:** Adequate ventilation and access to eye washing stations are required.

### Recommendations for Personal Protective Measures:

|                                  |   |
|----------------------------------|---|
| <b>Hand Protection:</b>          | Use gloves when handling this product. Inspect gloves prior to use for any visible damage. Use proper glove removal technique to prevent contact with skin. Dispose of gloves in accordance with applicable laws. Wash and dry hands after use. |
| <b>Eye Protection:</b>           | Use a face shield and/or safety goggles when handling this product. Use eyewear that has been approved by NIOSH.  |
| <b>Skin and Body Protection:</b> | Wear appropriate laboratory attire when handling this product, such as gloves, long pants, closed-toe shoes and a laboratory coat.  |

# Safety Data Sheet

Trade Name: DNA Wash Buffer

Revision Date: 4/29/2016

|                         |   |
|-------------------------|---|
| Respiratory Protection: | If engineering controls do not maintain airborne concentrations below recommended exposure limits (where applicable) or to an acceptable level (in countries where exposure limits have not been established), an approved respirator must be used. |
| Hygiene Measures:       | Avoid contact with eyes, skin and clothing. Wash hands thoroughly after using this product. Do not eat, drink or smoke when handling this product.  |

## Special Requirements for Personal Protective Equipment, Protective Clothing, Respirators, etc:

N/A

## Components with Workplace Control Parameters:

N/A

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES:

### Appearance:

|        |           |
|--------|-----------|
| Form:  | Liquid    |
| Color: | Colorless |

**Odor:** No data available

**Odor Threshold:** No data available

**pH:** 7.0-9.0

**Melting Point/Freezing Point:** No data available

**Initial Boiling Point/Boiling Range:** No data available

**Flash Point:** No data available

**Evaporation Rate:** No data available

**Flammability (Solid, Gas):** No data available

### Upper/Lower Limits on Flammability and Explosive Limits:

Flammability Limit (Upper): No data available

Flammability Limit (Lower): No data available

Explosive Limit (Upper): No data available

Explosive Limit (Lower): No data available

**Vapor Pressure:** No data available

**Vapor Density:** No data available

**Relative Density:** No data available.

# Safety Data Sheet

Trade Name: DNA Wash Buffer

Revision Date: 4/29/2016

**Solubility:**

Solubility in Water: No data available

Solubility (Other): No data available

**Partition Coefficient (n-octanol/  
water):**

No data available

**Auto-Ignition Temperature:**

No data available

**Decomposition Temperature:**

No data available

**Viscosity:**

No data available

## SECTION 10. STABILITY AND REACTIVITY:

**Reactivity:** No dangerous reactions known under conditions of normal use.

**Chemical Stability:** Stable under recommended storage conditions.

**Possibility of Hazardous Reactions:** N/A

**Conditions to Avoid:** N/A

**Materials to Avoid:** N/A

**Hazardous Decomposition Products:** N/A

## SECTION 11. TOXICOLOGICAL INFORMATION:

**Information on Likely Routes of Exposure:**

Inhalation: No data available

Ingestion: No data available

Skin: No data available

Eyes: No data available

**Information on Toxicological Effects:****Acute Toxicity:**

Oral LD50: No data available

Inhalation LC50: No data available

Dermal LD50: No data available

Other Information: No data available

Skin Corrosion/Irritation: No data available

Serious Eye Damage/Irritation: No data available

Respiratory or Skin Sensitization: No data available

# Safety Data Sheet

Trade Name: DNA Wash Buffer

Revision Date: 4/29/2016

Germ Cell Mutagenicity: No data available

## **Carcinogenicity:**

IARC: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH.

NTP: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.

OSHA: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.

Reproductive Toxicity: No data available

Teratogenicity: No data available

Specific Target Organ Toxicity - Single Exposure: No data available

Specific Target Organ Toxicity - Repeated Exposures: No data available

Aspiration Hazard: No data available

Signs and Symptoms of Exposure: No data available

Synergistic Effects: No data available

Additional Information: None known.

## **SECTION 12. ECOLOGICAL INFORMATION:**

Toxicity: No data available.

Persistence and Degradability: No data available.

Bioaccumulative Potential: No data available.

Mobility in Soil: No data available.

PBT and vPvB Assessment: No data available.

Other Adverse Effects: No data available.

# Safety Data Sheet

Trade Name: DNA Wash Buffer

Revision Date: 4/29/2016

## SECTION 13. DISPOSAL CONSIDERATIONS:

**Product:** Disposal should be in accordance with applicable regional, national and local laws and regulations.

Refer to Section 7: Handling and Storage and Section 8: Exposure Control/ Personal Protection for additional handling information and protection of employees.

**Contaminated Packaging:** Disposal must be made according to official regulations.

## SECTION 14. TRANSPORT INFORMATION:

**DOT:** Not regulated.

**IMDG:** Not regulated.

**IATA:** Not regulated.

## SECTION 15. REGULATORY INFORMATION:

### US Federal Regulations:

SARA 302 Components: None present or none present in regulated quantities.

SARA 313 Components: None present or none present in regulated quantities.

SARA 311/312 Hazards: No data available.

### US State Regulations:

California Prop. 65: No ingredient regulated by CA Prop 65 present.

Massachusetts RTK: No ingredient regulated by MS Right to Know Act present.

New Jersey & Community RTK: Water CAS No. 7732-18-5

Pennsylvania RTK: Water CAS No. 7732-18-5

Rhode Island RTK: No ingredient regulated by RI Right to Know Act present.



# Safety Data Sheet

Trade Name: DNA Wash Buffer

Revision Date: 4/29/2016

## SECTION 16. OTHER INFORMATION:

**Issue Date:** 9/10/2015

**Revision Date:** 4/29/2016

**Version:** 3.0

**Further Information:** No data available.

**Disclaimer:** The above information is believed to be correct but does not claim to be all inclusive. This document shall be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Omega Bio-tek, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

# Safety Data Sheet

Trade Name: Elution Buffer

Revision Date: 4/29/2016

## SECTION 1. IDENTIFICATION OF PRODUCT AND COMPANY:

**Product Name:** Elution Buffer

### **Manufacturer or Supplier Details:**

Company Name: Omega Bio-tek, Inc.  
Address: 400 Pinnacle Way, Suite 450  
Norcross, GA 30071  
USA  
  
Telephone: 1-800-832-8896  
Email: info@omegabiotek.com  
Emergency Telephone: CHEMTREC  
USA & Canada: 1-800-424-9300  
Outside USA & Canada: 1-703-527-3887

### **Recommended Use of the Chemical and Restrictions on Use:**

Recommended Use: For research use only.

## SECTION 2. HAZARDS IDENTIFICATION

### **Hazard Classification:**

Health Hazards: No known OSHA hazards.  
Not a dangerous substance according to GHS.

### **GHS Label Element:**

Hazard Symbol: N/A

### **Hazard Statements:**

Not a hazardous substance or mixture.

## SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS:

**Substance or Mixture:** Mixture of non-hazardous substances.

**Synonyms/Common Names:** N/A

# Safety Data Sheet

Trade Name: Elution Buffer

Revision Date: 4/29/2016

## SECTION 4. FIRST AID MEASURES:

|   |  |
|---|--|
| <b>General Advice:</b>  | Wear protective gloves/protective clothing/eye protection/face protection. |
| <b>If Inhaled:</b>  | If not breathing, give artificial respiration.                             |
| <b>In Case of Skin Contact:</b>   | Wash with soap and water.  |
| <b>In Case of Eye Contact:</b>  | Remove contact lenses, if present and easy to do. Continue rinsing.        |
| <b>If Swallowed:</b>  | N/A  |
| <b>Most Important Symptoms and Effects, Acute and Delayed:</b>                | N/A  |
| <b>Recommendations for Immediate Medical Attention and Special Treatment:</b> | Treat symptomatically. Symptoms may be delayed.                            |

## SECTION 5. FIREFIGHTING MEASURES:

|   |  |
|---|--|
| <b>Conditions of Flammability:</b>                    | Not flammable or combustible.  |
| <b>Suitable Extinguishing Media:</b>                  | Dry powder, CO <sub>2</sub> , water spray or regular foam.             |
| <b>Special Protective Equipment for Firefighters:</b> | Wear self-contained breathing apparatus, if necessary.                 |
| <b>Hazardous Combustion Products:</b>                 | N/A  |
| <b>Further Information:</b>                           | Prevent firefighting water from entering surface water or groundwater. |

# Safety Data Sheet

Trade Name: Elution Buffer

Revision Date: 4/29/2016

## SECTION 6. ACCIDENTAL RELEASE MEASURES:

|   |   |
|---|---|
| <b>Personal Precautions:</b>                                  | Avoid inhalation of vapors or mist. See Section 8: Exposure Controls/ Personal Protection for more information. |
| <b>Environmental Precautions:</b>                             | N/A   |
| <b>Emergency Procedures:</b>                                  | Evacuate personnel to safe areas and follow emergency response protocols.                                       |
| <b>Methods and Materials for Containment and Cleaning Up:</b> | Use inert absorbent material to soak up spill.  |

## SECTION 7. HANDLING AND STORAGE:

|                                       |  |
|---------------------------------------|--|
| <b>Precautions for Safe Handling:</b> | Handle in accordance with good chemical hygiene and safety practices. Protect from moisture. Avoid contact with eyes, skin and clothing. Avoid inhalation of vapors or mist. |
| <b>Conditions for Safe Storage:</b>   | Keep container tightly closed in a dry, well-ventilated place.   |
| <b>Incompatible Products:</b>         | N/A  |

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

**Appropriate Engineering Controls:** Adequate ventilation and access to eye washing stations are required.

### Recommendations for Personal Protective Measures:

|                                  |   |
|----------------------------------|---|
| <b>Hand Protection:</b>          | Use gloves when handling this product. Inspect gloves prior to use for any visible damage. Use proper glove removal technique to prevent contact with skin. Dispose of gloves in accordance with applicable laws. Wash and dry hands after use. |
| <b>Eye Protection:</b>           | Use a face shield and/or safety goggles when handling this product. Use eyewear that has been approved by NIOSH.  |
| <b>Skin and Body Protection:</b> | Wear appropriate laboratory attire when handling this product, such as gloves, long pants, closed-toe shoes and a laboratory coat.  |

# Safety Data Sheet

Trade Name: Elution Buffer

Revision Date: 4/29/2016

|                         |   |
|-------------------------|---|
| Respiratory Protection: | If engineering controls do not maintain airborne concentrations below recommended exposure limits (where applicable) or to an acceptable level (in countries where exposure limits have not been established), an approved respirator must be used. |
| Hygiene Measures:       | Avoid contact with eyes, skin and clothing. Wash hands thoroughly after using this product. Do not eat, drink or smoke when handling this product.  |

## Special Requirements for Personal Protective Equipment, Protective Clothing, Respirators, etc:

N/A

## Components with Workplace Control Parameters:

N/A

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES:

### Appearance:

|        |           |
|--------|-----------|
| Form:  | Liquid    |
| Color: | Colorless |

**Odor:** No data available

**Odor Threshold:** No data available

**pH:** 8.0-9.0

**Melting Point/Freezing Point:** No data available

**Initial Boiling Point/Boiling Range:** No data available

**Flash Point:** No data available

**Evaporation Rate:** No data available

**Flammability (Solid, Gas):** No data available

### Upper/Lower Limits on Flammability and Explosive Limits:

Flammability Limit (Upper): No data available

Flammability Limit (Lower): No data available

Explosive Limit (Upper): No data available

Explosive Limit (Lower): No data available

**Vapor Pressure:** No data available

**Vapor Density:** No data available

**Relative Density:** No data available.

# Safety Data Sheet

Trade Name: Elution Buffer

Revision Date: 4/29/2016

**Solubility:**

Solubility in Water: No data available

Solubility (Other): No data available

**Partition Coefficient (n-octanol/  
water):**

No data available

**Auto-Ignition Temperature:**

No data available

**Decomposition Temperature:**

No data available

**Viscosity:**

No data available

## SECTION 10. STABILITY AND REACTIVITY:

**Reactivity:** No dangerous reactions known under conditions of normal use.

**Chemical Stability:** Stable under recommended storage conditions.

**Possibility of Hazardous Reactions:** N/A

**Conditions to Avoid:** N/A

**Materials to Avoid:** N/A

**Hazardous Decomposition Products:** N/A

## SECTION 11. TOXICOLOGICAL INFORMATION:

**Information on Likely Routes of Exposure:**

Inhalation: No data available

Ingestion: No data available

Skin: No data available

Eyes: No data available

**Information on Toxicological Effects:****Acute Toxicity:**

Oral LD50: No data available

Inhalation LC50: No data available

Dermal LD50: No data available

Other Information: No data available

Skin Corrosion/Irritation: No data available

Serious Eye Damage/Irritation: No data available

Respiratory or Skin Sensitization: No data available

# Safety Data Sheet

Trade Name: Elution Buffer

Revision Date: 4/29/2016

Germ Cell Mutagenicity: No data available

## **Carcinogenicity:**

IARC: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH.

NTP: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.

OSHA: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.

Reproductive Toxicity: No data available

Teratogenicity: No data available

Specific Target Organ Toxicity - Single Exposure: No data available

Specific Target Organ Toxicity - Repeated Exposures: No data available

Aspiration Hazard: No data available

Signs and Symptoms of Exposure: No data available

Synergistic Effects: No data available

Additional Information: None known.

## **SECTION 12. ECOLOGICAL INFORMATION:**

Toxicity: No data available.

Persistence and Degradability: No data available.

Bioaccumulative Potential: No data available.

Mobility in Soil: No data available.

PBT and vPvB Assessment: No data available.

Other Adverse Effects: No data available.

# Safety Data Sheet

Trade Name: Elution Buffer

Revision Date: 4/29/2016

## SECTION 13. DISPOSAL CONSIDERATIONS:

**Product:** Disposal should be in accordance with applicable regional, national and local laws and regulations.

Refer to Section 7: Handling and Storage and Section 8: Exposure Control/ Personal Protection for additional handling information and protection of employees.

**Contaminated Packaging:** Disposal must be made according to official regulations.

## SECTION 14. TRANSPORT INFORMATION:

**DOT:** Not regulated.

**IMDG:** Not regulated.

**IATA:** Not regulated.

## SECTION 15. REGULATORY INFORMATION:

### US Federal Regulations:

SARA 302 Components: None present or none present in regulated quantities.

SARA 313 Components: None present or none present in regulated quantities.

SARA 311/312 Hazards: No data available.

### US State Regulations:

California Prop. 65: No ingredient regulated by CA Prop 65 present.

Massachusetts RTK: No ingredient regulated by MS Right to Know Act present.

New Jersey & Community RTK: Water CAS No. 7732-18-5

Pennsylvania RTK: Water CAS No. 7732-18-5

Rhode Island RTK: No ingredient regulated by RI Right to Know Act present.



# Safety Data Sheet

Trade Name: Elution Buffer

Revision Date: 4/29/2016

## SECTION 16. OTHER INFORMATION:

**Issue Date:** 9/10/2015

**Revision Date:** 4/29/2016

**Version:** 3.0

**Further Information:** No data available.

**Disclaimer:** The above information is believed to be correct but does not claim to be all inclusive. This document shall be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Omega Bio-tek, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

# Safety Data Sheet

Trade Name: HBC Buffer

Revision Date: 4/28/2016

## SECTION 1. IDENTIFICATION OF PRODUCT AND COMPANY:

**Product Name:** HBC Buffer

### Manufacturer or Supplier Details:

Company Name: Omega Bio-tek, Inc.  
Address: 400 Pinnacle Way, Suite 450  
Norcross, GA 30071  
USA  
  
Telephone: 1-800-832-8896  
Email: info@omegabiotek.com  
Emergency Telephone: CHEMTREC  
USA & Canada: 1-800-424-9300  
Outside USA & Canada: 1-703-527-3887

### Recommended Use of the Chemical and Restrictions on Use:

Recommended Use: For research use only.

## SECTION 2. HAZARDS IDENTIFICATION

### Hazard Classification:

|                 |   |   |
|-----------------|---|---|
| Health Hazards: | Acute toxicity (oral)<br>Skin corrosion/irritation<br>Serious eye damage/eye irritation<br>Specific target organ toxicity - single exposure | Category 4<br>Category 2<br>Category 2A<br>Category 3 |
|-----------------|---|---|

### GHS Label Element:

Hazard Symbol:



Signal Word:

Warning

### Hazard Statements:

H302 - Harmful if swallowed.  
H315 - Causes skin irritation.  
H319 - Causes serious skin irritation.  
H332 - Harmful if inhaled.

# Safety Data Sheet

Trade Name: HBC Buffer

Revision Date: 4/28/2016

## Precautionary Statements:

P261 - Avoid breathing dust/fumes/gas/mist/vapors/spray.  
P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P311 - Call a poison control center or doctor/physician.

## Potential Health Effects:

Aggravated Medical Conditions: May cause nervous system disturbances.  
Symptoms of Overexposure:  
Inhalation: Toxic if inhaled. May cause respiratory tract irritation.  
Skin: Harmful if absorbed through skin. May cause skin irritation.  
Eyes: Causes eye irritation.  
Ingestion: Harmful if swallowed.

## Carcinogenicity:

IARC No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.  
ACGIH No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH.  
OSHA No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.  
NTP No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.

## SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS:

Substance or Mixture: Mixture of the substances listed below with non-hazardous additions.

Synonyms/Common Names: Guanidinium chloride, Guanidinium monohydrochloride, Guanidinium hydrochloride

| Chemical Name  | CAS No. | Concentration (%) |
|--|---------|-------------------|
| Guanidine Hydrochloride<br>$\text{CH}_6\text{ClN}_3$ | 50-01-1 | 50-100%           |

# Safety Data Sheet

Trade Name: HBC Buffer

Revision Date: 4/28/2016

## SECTION 4. FIRST AID MEASURES:

|   |  |
|---|--|
| <b>General Advice:</b>  | Exit to a safe area. Consult a physician. Seek medical attention.                                      |
| <b>If Inhaled:</b>  | Move to fresh air. If not breathing, give artificial respiration. Seek medical attention.              |
| <b>In Case of Skin Contact:</b>   | Exit to safe area. Wash thoroughly with soap and water. Seek medical attention.                        |
| <b>In Case of Eye Contact:</b>  | Flush eyes with clean water for a minimum of 15 minutes, keeping eyes open. Seek medical attention.    |
| <b>If Swallowed:</b>  | Rinse mouth with water. Never give anything by mouth to an unconscious person. Seek medical attention. |
| <b>Most Important Symptoms and Effects, Acute and Delayed:</b>                | Irritating to eyes, respiratory system and skin. Narcotic effect.                                      |
| <b>Recommendations for Immediate Medical Attention and Special Treatment:</b> | Treat symptomatically. Symptoms may be delayed.  |

## SECTION 5. FIREFIGHTING MEASURES:

|   |  |
|---|--|
| <b>Conditions of Flammability:</b>                    | Not flammable or combustible.  |
| <b>Suitable Extinguishing Media:</b>                  | Dry powder, CO <sub>2</sub> , water spray or regular foam.                         |
| <b>Special Protective Equipment for Firefig ters:</b> | Wear self-contained breathing apparatus, if necessary.                             |
| <b>Hazardous Combustion Products:</b>                 | Carbon oxides (CO, CO <sub>2</sub> ), nitrogen oxides (NOx), hydrogen chloride gas |
| <b>Further Information:</b>                           | Prevent firefighting water from entering surface water or groundwater.             |

# Safety Data Sheet

Trade Name: HBC Buffer

Revision Date: 4/28/2016

## SECTION 6. ACCIDENTAL RELEASE MEASURES:

|   |   |
|---|---|
| <b>Personal Precautions:</b>                                  | Wear respiratory protection. Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. See Section 8: Exposure Controls/Personal Protection for more details. |
| <b>Environmental Precautions:</b>                             | Prevent further leakage if area is safe. Do not let product enter drains. Avoid discharge into the environment.   |
| <b>Emergency Procedures:</b>                                  | Evacuate personnel to safe areas and follow emergency response protocols.   |
| <b>Methods and Materials for Containment and Cleaning Up:</b> | Use inert absorbent material to soak up spill. Dispose of hazardous waste. Do not flush with water.   |

## SECTION 7. HANDLING AND STORAGE:

|                                       |  |
|---------------------------------------|--|
| <b>Precautions for Safe Handling:</b> | Handle in accordance with good chemical hygiene and safety practices. Protect from moisture. Avoid contact with eyes, skin and clothing. Avoid inhalation of vapors or mist. |
| <b>Conditions for Safe Storage:</b>   | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.  |
| <b>Incompatible Products:</b>         | Strong oxidizing agents  |

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

**Appropriate Engineering Controls:** Adequate ventilation and access to eye washing stations are required.

### Recommendations for Personal Protective Measures:

|                                  |   |
|----------------------------------|---|
| <b>Hand Protection:</b>          | Use gloves when handling this product. Inspect gloves prior to use for any visible damage. Use proper glove removal technique to prevent contact with skin. Dispose of gloves in accordance with applicable laws. Wash and dry hands after use. |
| <b>Eye Protection:</b>           | Use a face shield and/or safety goggles when handling this product. Use eyewear that has been approved by NIOSH.  |
| <b>Skin and Body Protection:</b> | Wear appropriate laboratory attire when handling this product, such as gloves, long pants, closed-toe shoes and a laboratory coat.  |

# Safety Data Sheet

Trade Name: HBC Buffer

Revision Date: 4/28/2016

|                         |   |
|-------------------------|---|
| Respiratory Protection: | If engineering controls do not maintain airborne concentrations below recommended exposure limits (where applicable) or to an acceptable level (in countries where exposure limits have not been established), an approved respirator must be used. |
| Hygiene Measures:       | Avoid contact with eyes, skin and clothing. Wash hands thoroughly after using this product. Do not eat, drink or smoke when handling this product.  |

**Special Requirements for Personal Protective Equipment, Protective Clothing, Respirators, etc:**

**Components with Workplace Control Parameters:**

| Components | CAS No. | Value | Control Parameters | Basis |
|------------|---------|-------|--------------------|-------|
|------------|---------|-------|--------------------|-------|

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES:

|   |                   |
|---|-------------------|
| <b>Appearance:</b>  |                   |
| Form:   | Liquid            |
| Color:  | Colorless         |
| Odor:   | No data available |
| Odor Threshold:   | No data available |
| pH:   | 3.0-5.0           |
| Melting Point/Freezing Point:                                   | No data available |
| Initial Boiling Point/Boiling Range:                            | No data available |
| Flash Point:  | No data available |
| Evaporation Rate:   | No data available |
| Flammability (Solid, Gas):                                      | No data available |
| <b>Upper/Lower Limits on Flammability and Explosive Limits:</b> |                   |
| Flammability Limit (Upper):                                     | No data available |
| Flammability Limit (Lower):                                     | No data available |
| Explosive Limit (Upper):  | No data available |
| Explosive Limit (Lower):  | No data available |
| Vapor Pressure:   | No data available |
| Vapor Density:  | No data available |

# Safety Data Sheet

Trade Name: HBC Buffer

Revision Date: 4/28/2016

**Relative Density:** No data available

**Solubility:**

Solubility in Water: No data available

Solubility (Other): No data available

**Partition Coefficient (n-octanol/water):** No data available

**Auto-Ignition Temperature:** No data available

**Decomposition Temperature:** No data available

**Viscosity:** No data available

## SECTION 10. STABILITY AND REACTIVITY:

**Reactivity:** No dangerous reaction known under conditions of normal use.

**Chemical Stability:** Stable under recommended storage conditions.

**Possibility of Hazardous Reactions:** Hazardous polymerization does not occur.

**Conditions to Avoid:** Contact with incompatible materials, excessive heat.

**Materials to Avoid:** Strong oxidizing agents, strong acids, inorganic salts, chlorinated compounds

**Hazardous Decomposition Products:** Carbon oxides, nitrogen oxides (NOx), hydrogen chloride gas

## SECTION 11. TOXICOLOGICAL INFORMATION:

### Information on Likely Routes of Exposure:

Inhalation: Harmful if inhaled. May cause respiratory tract irritation.

Ingestion: Harmful if swallowed.

Skin: Causes skin irritation.

Eyes: Causes serious eye irritation.

### Information on Toxicological Effects:

**Acute Toxicity:**

Oral LD50: 1,900 mg/kg (Rat)

Inhalation LC50: No data available

Dermal LD50: 2,000 mg/kg (Rabbit)

Other Information: No data available

# Safety Data Sheet

Trade Name: HBC Buffer

Revision Date: 4/28/2016

|                                    |  |
|------------------------------------|--|
| Skin Corrosion/Irritation:         | Causes skin irritation.  |
| Serious Eye Damage/Irritation:     | Causes serious eye irritation.   |
| Respiratory or Skin Sensitization: | Components of this product may be absorbed into the body through the skin. Irritating to skin and mucus membranes. |
| Germ Cell Mutagenicity:            | No mutagenic components identified.  |

## **Carcinogenicity:**

|        |  |
|--------|--|
| IARC:  | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.  |
| ACGIH: | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH. |
| NTP:   | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.   |
| OSHA:  | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.  |

**Reproductive Toxicity:** No component toxic to reproduction.

**Teratogenicity:** No data available.

**Specific Target Organ Toxicity - Single Exposure:** Respiratory tract irritation.

**Specific Target Organ Toxicity - Repeated Exposures:** No data available.

**Aspiration Hazard:** No data available.

**Signs and Symptoms of Exposure:** Irritating to eyes, respiratory system and skin. Narcotic effect.

**Synergistic Effects:** No data available.

**Additional Information:** None known.

## **SECTION 12. ECOLOGICAL INFORMATION:**

|                                       |  |
|---------------------------------------|--|
| <b>Toxicity:</b>                      | No data available.   |
| <b>Persistence and Degradability:</b> | No data available.   |
| <b>Bioaccumulative Potential:</b>     | No data available.   |
| <b>Mobility in Soil:</b>              | Product is water soluble and may spread in water systems.  |
| <b>PBT and vPvB Assessment:</b>       | No data available.   |
| <b>Other Adverse Effects:</b>         | The product components are not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment. |



# Safety Data Sheet

Trade Name: HBC Buffer

Revision Date: 4/28/2016

## SECTION 13. DISPOSAL CONSIDERATIONS:

**Product:**

Disposal should be in accordance with applicable regional, national and local laws and regulations.

S29 - Do not empty into drains.

S57 - Use appropriate container to avoid environmental contamination.

Refer to Section 7: Handling and Storage and Section 8: Exposure Control/Personal Protection for additional handling information and protection of employees.

**Contaminated Packaging:**

Disposal must be made according to official regulations.

To minimize exposure, please refer to Section 8: Exposure Control/Personal Protection of this SDS.

## SECTION 14. TRANSPORT INFORMATION:

**DOT:**

Not regulated.

**IMDG:**

Not regulated.

**IATA:**

Not regulated.

## SECTION 15. REGULATORY INFORMATION:

**US Federal Regulations:**

OSHA Hazards:

Aggravated Medical  
Conditions:

May cause nervous system disturbances.

Symptoms of  
Overexposure:

Inhalation:

Toxic if inhaled. May cause respiratory tract irritation.

Skin:

Harmful if absorbed through skin. May cause skin irritation.

Eyes:

Causes eye irritation.

Ingestion:

Harmful if swallowed.

SARA 302 Components:

None present or none present in regulated quantities.

SARA 313 Components:

None present or none present in regulated quantities.

SARA 311/312 Hazards:

Acute health hazard, chronic health hazard

# Safety Data Sheet

Trade Name: HBC Buffer

Revision Date: 4/28/2016

## US State Regulations:

|                             |  |
|-----------------------------|--|
| California Prop. 65:        | No ingredient regulated by CA Prop 65 present.           |
| New Jersey & Community RTK: | No ingredient regulated by NJ Right to Know Act present. |
| Massachusetts RTK:          | No ingredient regulated by MA Right to Know Act present. |
| Pennsylvania RTK:           | No ingredient regulated by PA Right to Know Act present. |
| Rhode Island RTK:           | No ingredient regulated by RI Right to Know Act present. |

## SECTION 16. OTHER INFORMATION:

|                      |  |
|----------------------|--|
| Issue Date:          | 8/19/2015  |
| Revision Date:       | 4/28/2016  |
| Version:             | 5.0  |
| Further Information: | No data available.   |
| Disclaimer:          | The above information is believed to be correct but does not claim to be all inclusive. This document shall be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Omega Bio-tek, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. |

# Safety Data Sheet

Trade Name: Lyticase

Revision Date: 6/12/2018

## SECTION 1. IDENTIFICATION OF PRODUCT AND COMPANY:

**Product Name:** Lyticase, from *Arthrobacter luteus*

### Manufacturer or Supplier Details:

**Company Name:** Bioland Scientific LLC  
**Address:** Bioland Scientific LLC  
14925 Paramount Blvd., Suite C  
Paramount, CA 90723  
USA  
**Telephone:** 1-562-377-2668  
**Fax:** 1-562-733-6008  
**Emergency Telephone:** 1-562-377-2668

### Recommended Use of the Chemical and Restrictions on Use:

**Recommended Use:** For research use only.

## SECTION 2. HAZARDS IDENTIFICATION

### Hazard Classification:

OSHA Hazards:

**Health Hazards:** Respiratory sensitization Category 1

### GHS Label Element:

**Hazard Symbol:**



**Signal Word:** Danger

**Hazard Statements:** H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

**Precautionary Statements:** P261 - Avoid breathing dust/fumes/gas/mist/vapors/spray.  
P285 - In case of inadequate ventilation, wear respiratory protection.  
P304 + P341 - IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.  
P342 + P311 - If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.  
P501 - Dispose of contents/container to an approved waste disposal plant.

# Safety Data Sheet

Trade Name: Lyticase

Revision Date: 6/12/2018

## SECTION 3. COMPOSITION OF/INFORMATION ON INGREDIENTS:

**Substance or Mixture:** Substances listed below.

**Synonyms/Common Names:** Zymolyase OEM equivalent

| Chemical Name | CAS No.    | Concentration (%) |
|---------------|------------|-------------------|
| Zymolyase     | 37340-57-1 | ≤ 100%            |

## SECTION 4. FIRST AID MEASURES:

**General Advice:** Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

**If Inhaled:** If breathed in, move person into fresh air. If not breathing, give artificial respiration.

**In Case of Skin Contact:** Wash off with soap and plenty of water. Consult a physician.

**In Case of Eye Contact:** Flush eyes with water as a precaution.

**If Swallowed:** Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

**Most Important Symptoms and Effects, Acute and Delayed:** The most important known symptoms and effects are described in the labelling (see Section 2.2) and/or in Section 11.

**Recommendations for Immediate Medical Attention and Special Treatment:** No data available

## SECTION 5. FIREFIGHTING MEASURES:

**Conditions of Flammability:** Not flammable or combustible.

**Suitable Extinguishing Media:** Water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

# Safety Data Sheet

Trade Name: Lyticase

Revision Date: 6/12/2018

**Special Protective Equipment for Firefighters:** Wear self-contained breathing apparatus, if necessary.

**Hazardous Combustion Products:** No data available

**Further Information:** Nature of decomposition products unknown.

## SECTION 6. ACCIDENTAL RELEASE MEASURES:

**Personal Precautions:** Use personal protective equipment. Avoid dust formation. Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust.

**Environmental Precautions:** Do not let product enter drains.

**Emergency Procedures:** Evacuate personnel to safe areas and follow emergency response protocols.

**Methods and Materials for Containment and Cleaning Up:** Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

## SECTION 7. HANDLING AND STORAGE:

**Precautions for Safe Handling:** Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Further processing of solid materials may result in the formation of combustible dusts. The potential for combustible dust formation should be taken into consideration before additional processing occurs.

Provide appropriate exhaust ventilation at places where dust is formed. For precautions, see Section 2.2.

**Conditions for Safe Storage:** Keep container tightly closed in a dry and well-ventilated place. Recommended storage temperature is -20°C. Storage class (TRGS 510): Non-combustible solids.

**Incompatible Products:** Strong oxidizing agents

# Safety Data Sheet

Trade Name: Lyticase

Revision Date: 6/12/2018

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

**Appropriate Engineering Controls:** Handle in accordance with good industrial hygiene and safety practices. Wash hands before breaks and at the end of the work day.

### Recommendations for Personal Protective Measures:

|                           |  |
|---------------------------|--|
| Hand Protection:          | Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.           |
| Eye Protection:           | Wear face shield and safety glasses. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU).  |
| Skin and Body Protection: | Complete suit protecting against chemicals. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.  |
| Respiratory Protection:   | For nuisance exposures, use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection, use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respiratory cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU). |
| Hygiene Measures:         | Avoid contact with skin, eyes and clothing. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.  |

### Special Requirements for Personal Protective Equipment, Protective Clothing, Respirators, etc:

### Components with Workplace Control Parameters:

Contains no substances with occupational exposure limit values.

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES:

### Appearance:

|        |                     |
|--------|---------------------|
| Form:  | Powder, lyophilized |
| Color: | Light brown         |

**Odor:** No data available

**Odor Threshold:** No data available

**pH:** No data available

**Melting Point/Freezing Point:** No data available

**Initial Boiling Point/Boiling Range:** No data available

# Safety Data Sheet

Trade Name: Lyticase

Revision Date: 6/12/2018

|   |                   |
|---|-------------------|
| <b>Flash Point:</b>   | No data available |
| <b>Evaporation Rate:</b>  | No data available |
| <b>Flammability (Solid, gas):</b>                               | No data available |
| <b>Upper/Lower Limits on Flammability and Explosive Limits:</b> |                   |
| Flammability Limit (Upper):                                     | No data available |
| Flammability Limit (Lower):                                     | No data available |
| Explosive Limit (Upper):  | No data available |
| Explosive Limit (Lower):  | No data available |
| <b>Vapor Pressure:</b>  | No data available |
| <b>Vapor Density:</b>   | No data available |
| <b>Relative Density:</b>  | No data available |
| <b>Solubility:</b>  |                   |
| Solubility in Water:  | No data available |
| Solubility (Other):   | No data available |
| <b>Partition Coefficient (n-octanol/<br/>water):</b>            | No data available |
| <b>Auto-Ignition Temperature:</b>                               | No data available |
| <b>Decomposition Temperature:</b>                               | No data available |
| <b>Viscosity:</b>   | No data available |

## SECTION 10. STABILITY AND REACTIVITY:

|  |  |
|--|--|
| <b>Reactivity:</b>                         | No data available                                |
| <b>Chemical Stability:</b>                 | Stable under recommended storage conditions.     |
| <b>Possibility of Hazardous Reactions:</b> | No data available                                |
| <b>Conditions to Avoid:</b>                | No data available                                |
| <b>Materials to Avoid:</b>                 | Strong oxidizing agents                          |
| <b>Hazardous Decomposition Products:</b>   | Other decomposition products - No data available |

## SECTION 11. TOXICOLOGICAL INFORMATION:

### Information on Likely Routes of Exposure:

|             |                   |
|-------------|-------------------|
| Inhalation: | No data available |
| Skin:       | No data available |
| Eyes:       | No data available |

# Safety Data Sheet

Trade Name: Lyticase

Revision Date: 6/12/2018

Ingestion: No data available

## Information on Toxicological Effects:

### Acute Toxicity:

Oral LD50: No data available

Inhalation LC50: No data available

Dermal LD50: No data available

Other Information: No data available

Skin Corrosion/Irritation: No data available

Serious Eye Damage/Irritation: No data available

Respiratory or Skin Sensitization: No data available

Germ Cell Mutagenicity: No data available

## Carcinogenicity:

IARC: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH.

NTP: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.

OSHA: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.

Reproductive Toxicity: No data available

Teratogenicity: No data available

Specific Target Organ Toxicity - Single Exposure: No data available

Specific Target Organ Toxicity - Repeated Exposures: No data available

Aspiration Hazard: No data available

Synergistic Effects: No data available

Additional Information: To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

## SECTION 12. ECOLOGICAL INFORMATION:

Toxicity: No data available.

Persistence and Degradability: No data available.

Bioaccumulative Potential: No data available.



# Safety Data Sheet

Trade Name: Lyticase

Revision Date: 6/12/2018

|                                 |   |
|---------------------------------|---|
| <b>Mobility in Soil:</b>        | No data available.  |
| <b>PBT and vPvB Assessment:</b> | PBT/vPvB assessment not available as chemical safety assessment not required/not conducted. |
| <b>Other Adverse Effects:</b>   | No data available   |

## SECTION 13. DISPOSAL CONSIDERATIONS:

|                                |  |
|--------------------------------|--|
| <b>Product:</b>                | Offer surplus and non-recyclable solutions to a licensed disposal company. Contact a licensed professional waste disposal service to dispose of this material. Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber. |
| <b>Contaminated Packaging:</b> | Dispose of as unused product.  |

## SECTION 14. TRANSPORT INFORMATION:

|              |                                 |
|--------------|---------------------------------|
| <b>DOT:</b>  | Not considered dangerous goods. |
| <b>IMDG:</b> | Not considered dangerous goods. |
| <b>IATA:</b> | Not considered dangerous goods. |

## SECTION 15. REGULATORY INFORMATION:

### US Federal Regulations:

|                       |   |
|-----------------------|---|
| SARA 302 Components:  | None present or none present in regulated quantities. |
| SARA 313 Components:  | None present or none present in regulated quantities. |
| SARA 311/312 Hazards: | Acute health hazard                                   |

### US State Regulations:

|                      |  |                |
|----------------------|--|----------------|
| California Prop. 65: | No ingredient regulated by CA Prop 65 present. |                |
| Massachusetts RTK:   | No ingredient regulated by MA RTK present.     |                |
| New Jersey RTK:      | <b>Component(s)</b>                            | <b>CAS-No.</b> |
|                      | Zymolyase                                      | 37340-57-1     |
| Pennsylvania RTK:    | Zymolyase                                      | 37340-57-1     |

# Safety Data Sheet

Trade Name: Lyticase

Revision Date: 6/12/2018

## SECTION 16. OTHER INFORMATION:

### Disclaimer:

\*\*\*\*\*

Bioland Scientific LLC provides the information contained herein in good faith but makes no representation as to its comprehensiveness or accuracy. This document is intended only as a guide to the appropriate precautionary handling of the material by a properly trained person using this product. Individuals receiving the information must exercise their independent judgment in determining its appropriateness for a particular purpose. BIOLAND SCIENTIFIC LLC MAKES NO REPRESENTATIONS OR WARRANTIES, EITHER EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE WITH RESPECT TO THE INFORMATION SET FORTH HEREIN OR THE PRODUCT TO WHICH THE INFORMATION REFERS. ACCORDINGLY, BIOLAND SCIENTIFIC LLC WILL NOT BE RESPONSIBLE FOR DAMAGES RESULTING FROM USE OF OR RELIANCE UPON THIS INFORMATION.

# Safety Data Sheet

Trade Name: Proteinase K Solution

Revision Date: 5/19/2016

## SECTION 1. IDENTIFICATION OF PRODUCT AND COMPANY:

**Product Name:** Proteinase K Solution

### Manufacturer or Supplier Details:

**Company Name:** Omega Bio-tek, Inc.  
**Address:** 400 Pinnacle Way, Suite 450  
Norcross, GA 30071  
USA  
**Telephone:** 1-800-832-8896  
**Email:** info@omegabiotek.com  
**Emergency Telephone:** CHEMTREC  
USA & Canada: 1-800-424-9300  
Outside USA & Canada: 1-703-527-3887

### Recommended Use of the Chemical and Restrictions on Use:

**Recommended Use:** For research use only.

## SECTION 2. HAZARDS IDENTIFICATION

### Hazard Classification:

|                 |                                   |            |
|-----------------|-----------------------------------|------------|
| OSHA Hazards:   | Respiratory sensitizer, irritant. |            |
| Health Hazards: | Skin irritation                   | Category 3 |
|                 | Respiratory sensitization         | Category 1 |

### GHS Label Element:

**Hazard Symbol:**



|                                  |  |
|----------------------------------|--|
| <b>Signal Word:</b>              | Danger   |
| <b>Hazard Statements:</b>        | H316 - Causes mild skin irritation.<br>H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.   |
| <b>Precautionary Statements:</b> | P261- Avoid breathing dust/fumes/gas/mist/vapors/spray.<br>P280 - Wear protective gloves/protective clothing/eye protection/face protection.<br>P284 - In case of inadequate ventilation, wear respiratory protection. |

# Safety Data Sheet

Trade Name: Proteinase K Solution

Revision Date: 5/19/2016

**Precautionary Statements:**

P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing.  
P332 + P 313 - If skin irritation occurs: Get medical advice/attention.  
P342 + P311 - If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.

**Potential Health Effects:**

Aggravated Medical Conditions:

Symptoms of Overexposure:

Inhalation:

May cause sensitization by inhalation.

Skin:

May cause mild skin irritation.

Eyes:

May cause eye irritation.

Ingestion:

May be harmful if swallowed. Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea.

**Carcinogenicity:**

IARC

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH.

OSHA

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.

NTP

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.

## SECTION 3. COMPOSITION OF/INFORMATION ON INGREDIENTS:

**Substance or Mixture:**

Mixture of the substances listed below with non-hazardous additions.

**Synonyms/Common Names:**

Endopeptidase K

| Chemical Name | CAS No.    | Concentration (%) |
|---------------|------------|-------------------|
| Proteinase K  | 39450-01-6 | 1-5%              |
| Glycerol      | 56-81-5    | 30-60%            |

# Safety Data Sheet

Trade Name: Proteinase K Solution

Revision Date: 5/19/2016

## SECTION 4. FIRST AID MEASURES:

|   |   |
|---|---|
| <b>General Advice:</b>  | Move out of dangerous area. Consult a physician. Show this safety data sheet to the doctor in attendance.       |
| <b>If Inhaled:</b>  | If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician. |
| <b>In Case of Skin Contact:</b>   | Wash off with soap and plenty of water. Consult a physician.  |
| <b>In Case of Eye Contact:</b>  | Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.                          |
| <b>If Swallowed:</b>  | Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.             |
| <b>Most Important Symptoms and Effects, Acute and Delayed:</b>                |   |
| <b>Recommendations for Immediate Medical Attention and Special Treatment:</b> | Treat symptomatically. Symptoms may be delayed.   |

## SECTION 5. FIREFIGHTING MEASURES:

|   |  |
|---|--|
| <b>Conditions of Flammability:</b>                    | Not flammable or combustible.  |
| <b>Suitable Extinguishing Media:</b>                  | Dry powder, CO <sub>2</sub> , water spray or alcohol-resistant foam. |
| <b>Special Protective Equipment for Firefighters:</b> | Wear self-contained breathing apparatus, if necessary.               |
| <b>Hazardous Combustion Products:</b>                 | N/A  |
| <b>Further Information:</b>                           |  |

# Safety Data Sheet

Trade Name: Proteinase K Solution

Revision Date: 5/19/2016

## SECTION 6. ACCIDENTAL RELEASE MEASURES:

|   |   |
|---|---|
| <b>Personal Precautions:</b>                                  | Use personal protective equipment. Avoid aerosolization. Avoid breathing vapors or mist. Ensure adequate ventilation. Evacuate personnel to safe areas.                               |
| <b>Environmental Precautions:</b>                             | Do not let product enter drains.  |
| <b>Emergency Procedures:</b>                                  | Evacuate personnel to safe areas and follow emergency response protocols.   |
| <b>Methods and Materials for Containment and Cleaning Up:</b> | Absorb on sand or vermiculite and place in closed containers for disposal. Pick up and transfer to properly labeled containers. Ensure proper ventilation and cleaning of spill site. |

## SECTION 7. HANDLING AND STORAGE:

|                                       |   |
|---------------------------------------|---|
| <b>Precautions for Safe Handling:</b> | Avoid contact with skin and eyes. Avoid inhalation of vapor or mist.                                      |
| <b>Conditions for Safe Storage:</b>   | Keep container tightly closed in a dry and well-ventilated place. Recommended storage temperature: 2-8°C. |
| <b>Incompatible Products:</b>         |   |

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

**Appropriate Engineering Controls:** Adequate ventilation and access to eye washing stations are required.

### Recommendations for Personal Protective Measures:

|                                  |   |
|----------------------------------|---|
| <b>Hand Protection:</b>          | Use gloves when handling this product. Inspect gloves prior to use for any visible damage. Use proper glove removal technique to prevent contact with skin. Dispose of gloves in accordance with applicable laws. Wash and dry hands after use. |
| <b>Eye Protection:</b>           | Use a face shield and/or safety goggles when handling this product. Use eyewear that has been approved by NIOSH.  |
| <b>Skin and Body Protection:</b> | Wear appropriate laboratory attire when handling this product, such as gloves, long pants, closed-toe shoes and a laboratory coat.  |

# Safety Data Sheet

Trade Name: Proteinase K Solution

Revision Date: 5/19/2016

Respiratory Protection: If engineering controls do not maintain airborne concentrations below recommended exposure limits (where applicable) or to an accepted level (in countries where exposure limits have not been established), an approved respirator must be used.

Hygiene Measures: Avoid contact with eyes, skin and clothing. Wash hands thoroughly after using this product. Do not eat, drink or smoke when handling this product.

## Special Requirements for Personal Protective Equipment, Protective Clothing, Respirators, etc:

### Components with Workplace Control Parameters:

| Components | CAS No. | Value                     | Control Parameters   | Basis   |
|------------|---------|---------------------------|----------------------|---|
| Glycerol   | 56-81-5 | PEL (Total dust)          | 15 mg/m <sup>3</sup> | USA. OSHA Table Z-1 Limits for Air Contaminants     |
|            |         | PEL (Respirable fraction) | 5 mg/m <sup>3</sup>  | USA. OSHA Table Z-1 Limits for Air Contaminants     |
|            |         | TWA                       | 5 mg/m <sup>3</sup>  | USA. Cal/OSHA Table Z-1 Limits for Air Contaminants |

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES:

### Appearance:

Form: Liquid

Color: Slightly yellow

Odor: No data available

Odor Threshold: No data available

pH: No data available

Melting Point/Freezing Point: No data available

Initial Boiling Point/Boiling Range: No data available

Flash Point: No data available

Evaporation Rate: No data available

Flammability (Solid, gas): No data available

# Safety Data Sheet

Trade Name: Proteinase K Solution

Revision Date: 5/19/2016

## Upper/Lower Limits on Flammability and Explosive Limits:

|                             |                   |
|-----------------------------|-------------------|
| Flammability Limit (Upper): | No data available |
| Flammability Limit (Lower): | No data available |
| Explosive Limit (Upper):    | No data available |
| Explosive Limit (Lower):    | No data available |

**Vapor Pressure:** No data available

**Vapor Density:** No data available

**Relative Density:** No data available

## Solubility:

|                      |                   |
|----------------------|-------------------|
| Solubility in Water: | No data available |
| Solubility (Other):  | No data available |

**Partition Coefficient (n-octanol/water):** No data available

**Auto-Ignition Temperature:** No data available

**Decomposition Temperature:** No data available

**Viscosity:** No data available

## SECTION 10. STABILITY AND REACTIVITY:

|  |   |
|--|---|
| <b>Reactivity:</b>                         | No dangerous reaction known under conditions of normal use. |
| <b>Chemical Stability:</b>                 | Stable under recommended storage conditions.                |
| <b>Possibility of Hazardous Reactions:</b> | N/A   |
| <b>Conditions to Avoid:</b>                | N/A   |
| <b>Materials to Avoid:</b>                 | N/A   |
| <b>Hazardous Decomposition Products:</b>   | N/A   |

## SECTION 11. TOXICOLOGICAL INFORMATION:

### Information on Likely Routes of Exposure:

|             |  |
|-------------|--|
| Inhalation: | May cause sensitization by inhalation.   |
| Skin:       | May cause mild skin irritation.  |
| Eyes:       | May cause eye irritation.  |
| Ingestion:  | May be harmful if swallowed. Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea. |



# Safety Data Sheet

Trade Name: Proteinase K Solution

Revision Date: 5/19/2016

## Information on Toxicological Effects:

| Acute Toxicity:                    | Proteinase K      | Glycerol                              |
|------------------------------------|-------------------|---------------------------------------|
| Oral LD50:                         | No data available | 12,600 mg/kg (rat)                    |
| Inhalation LC50:                   | No data available | > 570 mg/m <sup>3</sup> (rat)         |
| Dermal LD50:                       | No data available | > 10,000 mg/kg (rabbit)               |
| Other Information:                 | No data available |                                       |
| Skin Corrosion/Irritation:         | No data available | Mild skin irritation, 24 hrs (rabbit) |
| Serious Eye Damage/Irritation:     | No data available | Mild eye irritation, 24 hrs (rabbit)  |
| Respiratory or Skin Sensitization: | No data available |                                       |
| Germ Cell Mutagenicity:            | No data available |                                       |

## Carcinogenicity:

|        |  |
|--------|--|
| IARC:  | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.  |
| ACGIH: | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH. |
| NTP:   | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.   |
| OSHA:  | No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.  |

|  |                   |
|--|-------------------|
| Reproductive Toxicity:                               | No data available |
| Teratogenicity:                                      | No data available |
| Specific Target Organ Toxicity - Single Exposure:    | No data available |
| Specific Target Organ Toxicity - Repeated Exposures: | No data available |
| Aspiration Hazard:                                   | No data available |
| Signs and Symptoms of Exposure:                      |                   |
| Synergistic Effects:                                 | No data available |
| Additional Information:                              | No data available |

## SECTION 12. ECOLOGICAL INFORMATION:

|                        |   |
|------------------------|---|
| Toxicity:              | No data available.  |
| Toxicity to Fish LC50: | > 5,000 mg/L Carassius Auratus (goldfish)<br>44,000 mg/L Pimephales Promelas (fathead minnow)<br>67,500 mg/L Oncorhynchus Mykiss (rainbow trout)<br>> 100,100 mg/L Other fish |

# Safety Data Sheet

Trade Name: Proteinase K Solution

Revision Date: 5/19/2016

|                                       |  |
|---------------------------------------|--|
| <b>Persistence and Degradability:</b> | No data available.   |
| <b>Bioaccumulative Potential:</b>     | No data available.   |
| <b>Mobility in Soil:</b>              | Product is water soluble and may spread in water systems.  |
| <b>PBT and vPvB Assessment:</b>       | No data available.   |
| <b>Other Adverse Effects:</b>         | The product components are not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment. |

## SECTION 13. DISPOSAL CONSIDERATIONS:

|                                |   |
|--------------------------------|---|
| <b>Product:</b>                | Disposal should be in accordance with applicable regional, national and local laws and regulations.<br><br>Refer to Section 7: Handling and Storage and Section 8: Exposure Control/ Personal Protection for additional handling information and protection of employees. |
| <b>Contaminated Packaging:</b> | Disposal must be made according to official regulations.  |

## SECTION 14. TRANSPORT INFORMATION:

|              |                                 |
|--------------|---------------------------------|
| <b>DOT:</b>  | Not considered dangerous goods. |
| <b>IMDG:</b> | Not considered dangerous goods. |
| <b>IATA:</b> | Not considered dangerous goods. |

## SECTION 15. REGULATORY INFORMATION:

### US Federal Regulations:

|                       |   |
|-----------------------|---|
| OSHA Hazards:         | Respiratory sensitizer, irritant.                     |
| SARA 302 Components:  | None present or none present in regulated quantities. |
| SARA 313 Components:  | None present or none present in regulated quantities. |
| SARA 311/312 Hazards: | Chronic health hazard                                 |

### US State Regulations:

|                      |  |
|----------------------|--|
| California Prop. 65: | No ingredient regulated by CA Prop 65 present. |
| Massachusetts RTK:   | Glycerol 56-81-5                               |

# Safety Data Sheet

Trade Name: Proteinase K Solution

Revision Date: 5/19/2016

|                   |          |           |
|-------------------|----------|-----------|
| Pennsylvania RTK: | Water    | 7732-18-5 |
|                   | Glycerol | 56-81-5   |
| Rhode Island RTK: | Water    | 7732-18-5 |
|                   | Glycerol | 56-81-5   |

## SECTION 16. OTHER INFORMATION:

**Issue Date:** 9/8/2014

**Revision Date:** 5/19/2016

**Version:** 3.0

**Further Information:** No data available.

**Disclaimer:** The above information is believed to be correct but does not claim to be all inclusive. This document shall be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Omega Bio-tek, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

# Safety Data Sheet

Trade Name: RNase A

Revision Date: 12/30/2016

## SECTION 1. IDENTIFICATION OF PRODUCT AND COMPANY:

**Product Name:**

**Manufacturer or Supplier Details:**

|                      |  |
|----------------------|--|
| Company Name:        | Omega Bio-tek, Inc.  |
| Address:             | 400 Pinnacle Way, Suite 450<br>Norcross, GA 30071<br>USA                         |
| Telephone:           | 1-800-832-8896   |
| Email:               | info@omegabiotek.com   |
| Emergency Telephone: | CHEMTREC<br>USA & Canada: 1-800-424-9300<br>Outside USA & Canada: 1-703-527-3887 |

**Recommended Use of the Chemical and Restrictions on Use:**

|                  |                        |
|------------------|------------------------|
| Recommended Use: | For research use only. |
|------------------|------------------------|

## SECTION 2. HAZARDS IDENTIFICATION

**Hazard Classification:**

|                 |                        |            |
|-----------------|------------------------|------------|
| Health Hazards: | Respiratory sensitizer | Category 1 |
|-----------------|------------------------|------------|

**GHS Label Element:**



**Hazard Statements:**

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

**Precautionary Statements:**

P261 - Avoid breathing dust/fumes/gas/mist/spray.  
P280 - Wear protective gloves/protective clothing/eye protection/face protection.  
P284 - In case of inadequate ventilation, wear respiratory protection.  
P304 + P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.  
P332 + P313 - IF SKIN IRRITATION OCCURS: Get medical advice/attention.  
P342 + P311 - IF EXPERIENCING RESPIRATORY SYMPTOMS: Call a POISON CENTER or doctor/physician.

# Safety Data Sheet

Trade Name: RNase A

Revision Date: 12/30/2016

## SECTION 3. COMPOSITION OF/INFORMATION ON INGREDIENTS:

**Substance or Mixture:** Mixture of the substances listed below with non-hazardous additions.

**Synonyms/Common Names:** Ribonuclease A, ribonuclease I, ribonuclease 3' pyrimidinooligonucleotide hydrolase, pancreatic ribonuclease

| Chemical Name | CAS No.   | Concentration (%) |
|---------------|-----------|-------------------|
| Glycerol      | 56-81-5   | 30-60%            |
| RNase         | 9001-99-4 | 1-3%              |

## SECTION 4. FIRST AID MEASURES:

**General Advice:** Exit to a safe area. Consult a physician.

**If Inhaled:** If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

**In Case of Skin Contact:** Wash off with soap and plenty of water. Consult a physician.

**In Case of Eye Contact:** Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

**If Swallowed:** Rinse mouth with water. Never give anything by mouth to an unconscious person. Consult a physician.

**Most Important Symptoms and Effects, Acute and Delayed:** N/A

**Recommendations for Immediate Medical Attention and Special Treatment:** Treat symptomatically. Symptoms may be delayed.

## SECTION 5. FIREFIGHTING MEASURES:

**Conditions of Flammability:** Not flammable or combustible.

**Suitable Extinguishing Media:** Dry powder, CO<sub>2</sub>, water spray or alcohol-resistant foam.

# Safety Data Sheet

Trade Name: RNase A

Revision Date: 12/30/2016

**Special Protective Equipment for Firefighters:**

Wear self-contained breathing apparatus, if necessary.

**Hazardous Combustion Products:**

Carbon oxides, nitrogen oxides

**Further Information:**

Prevent firefighting water from entering surface water or groundwater.

## SECTION 6. ACCIDENTAL RELEASE MEASURES:

**Personal Precautions:**

Use personal protective equipment. Avoid breathing vapors or mist. Avoid contact with skin, eyes and clothing. Ensure adequate ventilation. See Section 8: Exposure Controls/Personal Protection for more information.

**Environmental Precautions:**

Do not let product enter drains.

**Emergency Procedures:**

Evacuate personnel to safe areas and follow emergency response protocols.

**Methods and Materials for Containment and Cleaning Up:**

Use inert material to soak up spill. Keep waste in suitable, closed containers for disposal.

## SECTION 7. HANDLING AND STORAGE:

**Precautions for Safe Handling:**

Handle in accordance with good chemical hygiene and safety practices. Avoid inhalation of vapor or mist.

**Conditions for Safe Storage:**

Store at 2-8°C. Keep containers tightly closed.

**Incompatible Products:**

Strong oxidizing agents

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

**Appropriate Engineering Controls:**

Adequate ventilation and access to eye washing stations are required.

# Safety Data Sheet

Trade Name: RNase A

Revision Date: 12/30/2016

## Recommendations for Personal Protective Measures:

|                           |   |
|---------------------------|---|
| Hand Protection:          | Use gloves when handling this product. Inspect gloves prior to use for any visible damage. Use proper glove removal technique to prevent contact with skin. Dispose of gloves in accordance with applicable laws. Wash and dry hands after use.     |
| Eye Protection:           | Use a face shield and/or safety goggles when handling this product. Use eyewear that has been approved by NIOSH.  |
| Skin and Body Protection: | Wear appropriate laboratory attire when handling this product, such as gloves, long pants, closed-toe shoes and a laboratory coat.  |
| Respiratory Protection:   | If engineering controls do not maintain airborne concentrations below recommended exposure limits (where applicable) or to an acceptable level (in countries where exposure limits have not been established), an approved respirator must be used. |
| Hygiene Measures:         | Avoid contact with eyes, skin and clothing. Wash hands thoroughly after using this product. Do not eat, drink or smoke when handling this product.  |

## Special Requirements for Personal Protective Equipment, Protective Clothing, Respirators, etc:

N/A

## Components with Workplace Control Parameters:

| Components | CAS No. | Value | Control Parameters   | Basis   |
|------------|---------|-------|----------------------|---|
| Glycerol   | 56-81-5 | PEL   | 15 mg/m <sup>3</sup> | USA. OSHA Table Z-1 Limits for Air Contaminants |
|            |         | PEL   | 5 mg/m <sup>3</sup>  | USA. OSHA Table Z-1 Limits for Air Contaminants |
|            |         | TWA   | 10 mg/m <sup>3</sup> | USA. ACGIH Threshold Limit Values (TLV)         |

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES:

### Appearance:

|                                      |                   |
|--------------------------------------|-------------------|
| Form:                                | Liquid            |
| Color:                               | Colorless         |
| Odor:                                | No data available |
| Odor Threshold:                      | No data available |
| pH:                                  | No data available |
| Melting Point/Freezing Point:        | No data available |
| Initial Boiling Point/Boiling Range: | No data available |
| Flash Point:                         | No data available |

# Safety Data Sheet

Trade Name: RNase A

Revision Date: 12/30/2016

|   |                   |
|---|-------------------|
| <b>Evaporation Rate:</b>  | No data available |
| <b>Flammability (Solid, Gas):</b>                               | No data available |
| <b>Upper/Lower Limits on Flammability and Explosive Limits:</b> |                   |
| Flammability Limit (Upper):                                     | No data available |
| Flammability Limit (Lower):                                     | No data available |
| <b>Vapor Pressure:</b>  | No data available |
| <b>Vapor Density:</b>   | No data available |
| <b>Relative Density:</b>  | No data available |
| <b>Solubility:</b>  |                   |
| Solubility in Water:  | No data available |
| Solubility (Other):   | No data available |
| <b>Partition Coefficient (n-octanol/<br/>water):</b>            | No data available |
| <b>Auto-Ignition Temperature:</b>                               | No data available |
| <b>Decomposition Temperature:</b>                               | No data available |
| <b>Viscosity:</b>   | No data available |

## SECTION 10. STABILITY AND REACTIVITY:

|  |  |
|--|--|
| <b>Reactivity:</b>                         | No dangerous reaction known under conditions of normal use |
| <b>Chemical Stability:</b>                 | Stable under recommended storage conditions                |
| <b>Possibility of Hazardous Reactions:</b> | None under normal processing                               |
| <b>Conditions to Avoid:</b>                | Extreme temperatures and direct sunlight                   |
| <b>Materials to Avoid:</b>                 | Strong oxidizing agents                                    |
| <b>Hazardous Decomposition Products:</b>   | No data available  |

## SECTION 11. TOXICOLOGICAL INFORMATION:

### Information on Toxicological Effects:

| Acute Toxicity:                | RNase A           | Glycerol                              |
|--------------------------------|-------------------|---------------------------------------|
| Oral LD50:                     | No data available | 12,600 mg/kg (rat)                    |
| Inhalation LC50:               | No data available | No data available                     |
| Dermal LD50:                   | No data available | > 10,000 mg/kg (rabbit)               |
| Other Information:             | No data available | No data available                     |
| Skin Corrosion/Irritation:     | No data available | Mild skin irritation, 24 hrs (rabbit) |
| Serious Eye Damage/Irritation: | No data available | Mild eye irritation, 24 hrs (rabbit)  |



# Safety Data Sheet

Trade Name: RNase A

Revision Date: 12/30/2016

Respiratory or Skin Sensitization: No data available

Germ Cell Mutagenicity: No data available

## **Carcinogenicity:**

IARC: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH.

NTP: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.

OSHA: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.

**Reproductive Toxicity:** No data available

**Teratogenicity:** No data available

**Specific Target Organ Toxicity - Single Exposure:** No data available

**Specific Target Organ Toxicity - Repeated Exposures:** No data available

**Aspiration Hazard:** No data available

**Signs and Symptoms of Exposure:** Prolonged or repeated exposure may cause headache, nausea, vomiting.

**Synergistic Effects:** No data available

**Additional Information:** Kidney: Irregularities - Based on human evidence (Glycerol)

## **SECTION 12. ECOLOGICAL INFORMATION:**

**Toxicity:** No data available.

**Persistence and Degradability:** No data available.

**Bioaccumulative Potential:** No data available.

**Mobility in Soil:** No data available.

**PBT and vPvB Assessment:** No data available.

**Other Adverse Effects:** No data available.

# Safety Data Sheet

Trade Name: RNase A

Revision Date: 12/30/2016

## SECTION 13. DISPOSAL CONSIDERATIONS:

**Product:**

Disposal should be in accordance with applicable regional, national and local laws and regulations.

S29 - Do not empty into drains.

S57 - Use appropriate container to avoid environmental contamination.

Refer to Section 7: Handling and Storage and Section 8: Exposure Control/Personal Protection for additional handling information and protection of employees.

**Contaminated Packaging:**

Disposal must be made according to official regulations.

To minimize exposure, please refer to Section 8: Exposure Control/Personal Protection of this SDS.

## SECTION 14. TRANSPORT INFORMATION:

**DOT:**

Not considered dangerous goods.

**IMDG:**

Not considered dangerous goods.

**IATA:**

Not considered dangerous goods.

## SECTION 15. REGULATORY INFORMATION:

**US Federal Regulations:**

SARA 302 Components:

None present or none present in regulated quantities.

SARA 313 Components:

None present or none present in regulated quantities.

SARA 311/312 Hazards:

Acute health hazard

**US State Regulations:**

California Prop. 65:

No ingredient regulated by CA Prop 65 present.

# Safety Data Sheet

Trade Name: RNase A

Revision Date: 12/30/2016

## SECTION 16. OTHER INFORMATION:

**Issue Date:** 9/9/2014

**Revision Date:** 12/30/2016

**Version:** 2.5

**Further Information:** No data available.

**Disclaimer:** The above information is believed to be correct but does not claim to be all inclusive. This document shall be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Omega Bio-tek, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

# Safety Data Sheet

Trade Name: SE Buffer

Revision Date: 6/29/2018

## SECTION 1. IDENTIFICATION OF PRODUCT AND COMPANY:

**Product Name:** SE Buffer

### **Manufacturer or Supplier Details:**

|                      |  |
|----------------------|--|
| Company Name:        | Bioland Scientific LLC   |
| Address:             | Bioland Scientific LLC<br>14925 Paramount Blvd., Suite C<br>Paramount, CA 90723<br>USA |
| Telephone:           | 1-562-377-2668   |
| Email:               | 1-562-733-6008   |
| Emergency Telephone: | 1-562-377-2668   |

### **Recommended Use of the Chemical and Restrictions on Use:**

|                  |                        |
|------------------|------------------------|
| Recommended Use: | For research use only. |
|------------------|------------------------|

## SECTION 2. HAZARDS IDENTIFICATION

### **Hazard Classification:**

|                 |   |
|-----------------|---|
| Health Hazards: | No known OSHA hazards.<br>Not a dangerous substance according to GHS. |
|-----------------|---|

### **GHS Label Element:**

|                |     |
|----------------|-----|
| Hazard Symbol: | N/A |
|----------------|-----|

|                           |                                       |
|---------------------------|---------------------------------------|
| <b>Hazard Statements:</b> | Not a hazardous substance or mixture. |
|---------------------------|---------------------------------------|

## SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS:

|                              |                                      |
|------------------------------|--------------------------------------|
| <b>Substance or Mixture:</b> | Mixture of non-hazardous substances. |
|------------------------------|--------------------------------------|

|                               |     |
|-------------------------------|-----|
| <b>Synonyms/Common Names:</b> | N/A |
|-------------------------------|-----|

# Safety Data Sheet

Trade Name: SE Buffer

Revision Date: 6/29/2018

## SECTION 4. FIRST AID MEASURES:

|   |  |
|---|--|
| <b>General Advice:</b>  | Wear protective gloves/protective clothing/eye protection/face protection. |
| <b>If Inhaled:</b>  | If not breathing, give artificial respiration.                             |
| <b>In Case of Skin Contact:</b>   | Wash with soap and water.  |
| <b>In Case of Eye Contact:</b>  | Remove contact lenses, if present and easy to do. Continue rinsing.        |
| <b>If Swallowed:</b>  | N/A  |
| <b>Most Important Symptoms and Effects, Acute and Delayed:</b>                | N/A  |
| <b>Recommendations for Immediate Medical Attention and Special Treatment:</b> | Treat symptomatically. Symptoms may be delayed.                            |

## SECTION 5. FIREFIGHTING MEASURES:

|   |  |
|---|--|
| <b>Conditions of Flammability:</b>                    | Not flammable or combustible.  |
| <b>Suitable Extinguishing Media:</b>                  | Dry powder, CO <sub>2</sub> , water spray or regular foam.             |
| <b>Special Protective Equipment for Firefighters:</b> | Wear self-contained breathing apparatus, if necessary.                 |
| <b>Hazardous Combustion Products:</b>                 | N/A  |
| <b>Further Information:</b>                           | Prevent firefighting water from entering surface water or groundwater. |

# Safety Data Sheet

Trade Name: SE Buffer

Revision Date: 6/29/2018

## SECTION 6. ACCIDENTAL RELEASE MEASURES:

|   |   |
|---|---|
| <b>Personal Precautions:</b>                                  | Avoid inhalation of vapors or mist. See Section 8: Exposure Controls/ Personal Protection for more information. |
| <b>Environmental Precautions:</b>                             | N/A   |
| <b>Emergency Procedures:</b>                                  | Evacuate personnel to safe areas and follow emergency response protocols.                                       |
| <b>Methods and Materials for Containment and Cleaning Up:</b> | Use inert absorbent material to soak up spill.  |

## SECTION 7. HANDLING AND STORAGE:

|                                       |  |
|---------------------------------------|--|
| <b>Precautions for Safe Handling:</b> | Handle in accordance with good chemical hygiene and safety practices. Protect from moisture. Avoid contact with eyes, skin and clothing. Avoid inhalation of vapors or mist. |
| <b>Conditions for Safe Storage:</b>   | Keep container tightly closed in a dry, well-ventilated place.   |
| <b>Incompatible Products:</b>         | N/A  |

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

**Appropriate Engineering Controls:** Adequate ventilation and access to eye washing stations are required.

### Recommendations for Personal Protective Measures:

|                                  |   |
|----------------------------------|---|
| <b>Hand Protection:</b>          | Use gloves when handling this product. Inspect gloves prior to use for any visible damage. Use proper glove removal technique to prevent contact with skin. Dispose of gloves in accordance with applicable laws. Wash and dry hands after use. |
| <b>Eye Protection:</b>           | Use a face shield and/or safety goggles when handling this product. Use eyewear that has been approved by NIOSH.  |
| <b>Skin and Body Protection:</b> | Wear appropriate laboratory attire when handling this product, such as gloves, long pants, closed-toe shoes and a laboratory coat.  |

# Safety Data Sheet

Trade Name: SE Buffer

Revision Date: 6/29/2018

|                         |   |
|-------------------------|---|
| Respiratory Protection: | If engineering controls do not maintain airborne concentrations below recommended exposure limits (where applicable) or to an acceptable level (in countries where exposure limits have not been established), an approved respirator must be used. |
| Hygiene Measures:       | Avoid contact with eyes, skin and clothing. Wash hands thoroughly after using this product. Do not eat, drink or smoke when handling this product.  |

## Special Requirements for Personal Protective Equipment, Protective Clothing, Respirators, etc:

N/A

## Components with Workplace Control Parameters:

N/A

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES:

### Appearance:

|        |           |
|--------|-----------|
| Form:  | Liquid    |
| Color: | Colorless |

**Odor:** No data available

**Odor Threshold:** No data available

**pH:** No data available

**Melting Point/Freezing Point:** No data available

**Initial Boiling Point/Boiling Range:** No data available

**Flash Point:** No data available

**Evaporation Rate:** No data available

**Flammability (Solid, Gas):** No data available

### Upper/Lower Limits on Flammability and Explosive Limits:

Flammability Limit (Upper): No data available

Flammability Limit (Lower): No data available

Explosive Limit (Upper): No data available

Explosive Limit (Lower): No data available

**Vapor Pressure:** No data available

**Vapor Density:** No data available

**Relative Density:** No data available.

# Safety Data Sheet

Trade Name: SE Buffer

Revision Date: 6/29/2018

**Solubility:**

Solubility in Water: No data available

Solubility (Other): No data available

**Partition Coefficient (n-octanol/  
water):**

No data available

**Auto-Ignition Temperature:**

No data available

**Decomposition Temperature:**

No data available

**Viscosity:**

No data available

## SECTION 10. STABILITY AND REACTIVITY:

**Reactivity:** No dangerous reactions known under conditions of normal use.

**Chemical Stability:** Stable under recommended storage conditions.

**Possibility of Hazardous Reactions:** N/A

**Conditions to Avoid:** N/A

**Materials to Avoid:** N/A

**Hazardous Decomposition Products:** N/A

## SECTION 11. TOXICOLOGICAL INFORMATION:

**Information on Likely Routes of Exposure:**

Inhalation: No data available

Ingestion: No data available

Skin: No data available

Eyes: No data available

**Information on Toxicological Effects:****Acute Toxicity:**

Oral LD50: No data available

Inhalation LC50: No data available

Dermal LD50: No data available

Other Information: No data available

Skin Corrosion/Irritation: No data available

Serious Eye Damage/Irritation: No data available

Respiratory or Skin Sensitization: No data available



# Safety Data Sheet

Trade Name: SE Buffer

Revision Date: 6/29/2018

Germ Cell Mutagenicity: No data available

## **Carcinogenicity:**

IARC: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH.

NTP: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.

OSHA: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.

Reproductive Toxicity: No data available

Teratogenicity: No data available

Specific Target Organ Toxicity - Single Exposure: No data available

Specific Target Organ Toxicity - Repeated Exposures: No data available

Aspiration Hazard: No data available

Signs and Symptoms of Exposure: No data available

Synergistic Effects: No data available

Additional Information: None known.

## **SECTION 12. ECOLOGICAL INFORMATION:**

Toxicity: No data available.

Persistence and Degradability: No data available.

Bioaccumulative Potential: No data available.

Mobility in Soil: No data available.

PBT and vPvB Assessment: No data available.

Other Adverse Effects: No data available.

# Safety Data Sheet

Trade Name: SE Buffer

Revision Date: 6/29/2018

## SECTION 13. DISPOSAL CONSIDERATIONS:

**Product:** Disposal should be in accordance with applicable regional, national and local laws and regulations.

Refer to Section 7: Handling and Storage and Section 8: Exposure Control/ Personal Protection for additional handling information and protection of employees.

**Contaminated Packaging:** Disposal must be made according to official regulations.

## SECTION 14. TRANSPORT INFORMATION:

**DOT:** Not regulated.

**IMDG:** Not regulated.

**IATA:** Not regulated.

## SECTION 15. REGULATORY INFORMATION:

### US Federal Regulations:

SARA 302 Components: None present or none present in regulated quantities.

SARA 313 Components: None present or none present in regulated quantities.

SARA 311/312 Hazards: No data available.

### US State Regulations:

California Prop. 65: No ingredient regulated by CA Prop 65 present.

Massachusetts RTK: No ingredient regulated by MS Right to Know Act present.

New Jersey & Community RTK: Water CAS No. 7732-18-5

Pennsylvania RTK: Water CAS No. 7732-18-5

Rhode Island RTK: No ingredient regulated by RI Right to Know Act present.

# Safety Data Sheet

Trade Name: SE Buffer

Revision Date: 6/29/2018

## SECTION 16. OTHER INFORMATION:

### Disclaimer:

\*\*\*\*\*

Bioland Scientific LLC provides the information contained herein in good faith but makes no representation as to its comprehensiveness or accuracy. This document is intended only as a guide to the appropriate precautionary handling of the material by a properly trained person using this product. Individuals receiving the information must exercise their independent judgment in determining its appropriateness for a particular purpose. BIOLAND SCIENTIFIC LLC MAKES NO REPRESENTATIONS OR WARRANTIES, EITHER EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE WITH RESPECT TO THE INFORMATION SET FORTH HEREIN OR THE PRODUCT TO WHICH THE INFORMATION REFERS. ACCORDINGLY, BIOLAND SCIENTIFIC LLC WILL NOT BE RESPONSIBLE FOR DAMAGES RESULTING FROM USE OF OR RELIANCE UPON THIS INFORMATION.

# Safety Data Sheet

Trade Name: TL Buffer

Revision Date: 4/29/2016

## SECTION 1. IDENTIFICATION OF PRODUCT AND COMPANY:

**Product Name:** TL Buffer

### **Manufacturer or Supplier Details:**

**Company Name:** Omega Bio-tek, Inc.  
**Address:** 400 Pinnacle Way, Suite 450  
Norcross, GA 30071  
USA  
**Telephone:** 1-800-832-8896  
**Email:** info@omegabiotek.com  
**Emergency Telephone:** CHEMTREC  
USA & Canada: 1-800-424-9300  
Outside USA & Canada: 1-703-527-3887

### **Recommended Use of the Chemical and Restrictions on Use:**

**Recommended Use:** For research use only.

## SECTION 2. HAZARDS IDENTIFICATION

### **Hazard Classification:**

**Health Hazards:** No known OSHA hazards.  
Not a dangerous substance according to GHS.

### **GHS Label Element:**

**Hazard Symbol:** N/A

**Hazard Statements:** Not a hazardous substance or mixture.

## SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS:

**Substance or Mixture:** Mixture of non-hazardous substances.

**Synonyms/Common Names:** N/A

# Safety Data Sheet

Trade Name: TL Buffer

Revision Date: 4/29/2016

## SECTION 4. FIRST AID MEASURES:

|   |  |
|---|--|
| <b>General Advice:</b>  | Wear protective gloves/protective clothing/eye protection/face protection. |
| <b>If Inhaled:</b>  | If not breathing, give artificial respiration.                             |
| <b>In Case of Skin Contact:</b>   | Wash with soap and water.  |
| <b>In Case of Eye Contact:</b>  | Remove contact lenses, if present and easy to do. Continue rinsing.        |
| <b>If Swallowed:</b>  | N/A  |
| <b>Most Important Symptoms and Effects, Acute and Delayed:</b>                | N/A  |
| <b>Recommendations for Immediate Medical Attention and Special Treatment:</b> | Treat symptomatically. Symptoms may be delayed.                            |

## SECTION 5. FIREFIGHTING MEASURES:

|   |  |
|---|--|
| <b>Conditions of Flammability:</b>                    | Not flammable or combustible.  |
| <b>Suitable Extinguishing Media:</b>                  | Dry powder, CO <sub>2</sub> , water spray or regular foam.             |
| <b>Special Protective Equipment for Firefighters:</b> | Wear self-contained breathing apparatus, if necessary.                 |
| <b>Hazardous Combustion Products:</b>                 | N/A  |
| <b>Further Information:</b>                           | Prevent firefighting water from entering surface water or groundwater. |

# Safety Data Sheet

Trade Name: TL Buffer

Revision Date: 4/29/2016

## SECTION 6. ACCIDENTAL RELEASE MEASURES:

|   |   |
|---|---|
| <b>Personal Precautions:</b>                                  | Avoid inhalation of vapors or mist. See Section 8: Exposure Controls/ Personal Protection for more information. |
| <b>Environmental Precautions:</b>                             | N/A   |
| <b>Emergency Procedures:</b>                                  | Evacuate personnel to safe areas and follow emergency response protocols.                                       |
| <b>Methods and Materials for Containment and Cleaning Up:</b> | Use inert absorbent material to soak up spill.  |

## SECTION 7. HANDLING AND STORAGE:

|                                       |  |
|---------------------------------------|--|
| <b>Precautions for Safe Handling:</b> | Handle in accordance with good chemical hygiene and safety practices. Protect from moisture. Avoid contact with eyes, skin and clothing. Avoid inhalation of vapors or mist. |
| <b>Conditions for Safe Storage:</b>   | Keep container tightly closed in a dry, well-ventilated place.   |
| <b>Incompatible Products:</b>         | N/A  |

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

**Appropriate Engineering Controls:** Adequate ventilation and access to eye washing stations are required.

### Recommendations for Personal Protective Measures:

|                                  |   |
|----------------------------------|---|
| <b>Hand Protection:</b>          | Use gloves when handling this product. Inspect gloves prior to use for any visible damage. Use proper glove removal technique to prevent contact with skin. Dispose of gloves in accordance with applicable laws. Wash and dry hands after use. |
| <b>Eye Protection:</b>           | Use a face shield and/or safety goggles when handling this product. Use eyewear that has been approved by NIOSH.  |
| <b>Skin and Body Protection:</b> | Wear appropriate laboratory attire when handling this product, such as gloves, long pants, closed-toe shoes and a laboratory coat.  |

# Safety Data Sheet

Trade Name: TL Buffer

Revision Date: 4/29/2016

|                         |   |
|-------------------------|---|
| Respiratory Protection: | If engineering controls do not maintain airborne concentrations below recommended exposure limits (where applicable) or to an acceptable level (in countries where exposure limits have not been established), an approved respirator must be used. |
| Hygiene Measures:       | Avoid contact with eyes, skin and clothing. Wash hands thoroughly after using this product. Do not eat, drink or smoke when handling this product.  |

## Special Requirements for Personal Protective Equipment, Protective Clothing, Respirators, etc:

N/A

## Components with Workplace Control Parameters:

N/A

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES:

### Appearance:

|        |           |
|--------|-----------|
| Form:  | Liquid    |
| Color: | Colorless |

**Odor:** No data available

**Odor Threshold:** No data available

**pH:** 7.0-9.0

**Melting Point/Freezing Point:** No data available

**Initial Boiling Point/Boiling Range:** No data available

**Flash Point:** No data available

**Evaporation Rate:** No data available

**Flammability (Solid, Gas):** No data available

### Upper/Lower Limits on Flammability and Explosive Limits:

Flammability Limit (Upper): No data available

Flammability Limit (Lower): No data available

Explosive Limit (Upper): No data available

Explosive Limit (Lower): No data available

**Vapor Pressure:** No data available

**Vapor Density:** No data available

**Relative Density:** No data available.

# Safety Data Sheet

Trade Name: TL Buffer

Revision Date: 4/29/2016

**Solubility:**

Solubility in Water: No data available

Solubility (Other): No data available

**Partition Coefficient (n-octanol/  
water):**

No data available

**Auto-Ignition Temperature:**

No data available

**Decomposition Temperature:**

No data available

**Viscosity:**

No data available

## SECTION 10. STABILITY AND REACTIVITY:

**Reactivity:** No dangerous reactions known under conditions of normal use.

**Chemical Stability:** Stable under recommended storage conditions.

**Possibility of Hazardous Reactions:** N/A

**Conditions to Avoid:** N/A

**Materials to Avoid:** N/A

**Hazardous Decomposition Products:** N/A

## SECTION 11. TOXICOLOGICAL INFORMATION:

**Information on Likely Routes of Exposure:**

Inhalation: No data available

Ingestion: No data available

Skin: No data available

Eyes: No data available

**Information on Toxicological Effects:****Acute Toxicity:**

Oral LD50: No data available

Inhalation LC50: No data available

Dermal LD50: No data available

Other Information: No data available

Skin Corrosion/Irritation: No data available

Serious Eye Damage/Irritation: No data available

Respiratory or Skin Sensitization: No data available



# Safety Data Sheet

Trade Name: TL Buffer

Revision Date: 4/29/2016

Germ Cell Mutagenicity: No data available

## **Carcinogenicity:**

IARC: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by ACGIH.

NTP: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by NTP.

OSHA: No component of this product that is present at 0.1% is identified as probable, possible or confirmed human carcinogen by OSHA.

Reproductive Toxicity: No data available

Teratogenicity: No data available

Specific Target Organ Toxicity - Single Exposure: No data available

Specific Target Organ Toxicity - Repeated Exposures: No data available

Aspiration Hazard: No data available

Signs and Symptoms of Exposure: No data available

Synergistic Effects: No data available

Additional Information: None known.

## **SECTION 12. ECOLOGICAL INFORMATION:**

Toxicity: No data available.

Persistence and Degradability: No data available.

Bioaccumulative Potential: No data available.

Mobility in Soil: No data available.

PBT and vPvB Assessment: No data available.

Other Adverse Effects: No data available.

# Safety Data Sheet

Trade Name: TL Buffer

Revision Date: 4/29/2016

## SECTION 13. DISPOSAL CONSIDERATIONS:

**Product:** Disposal should be in accordance with applicable regional, national and local laws and regulations.

Refer to Section 7: Handling and Storage and Section 8: Exposure Control/ Personal Protection for additional handling information and protection of employees.

**Contaminated Packaging:** Disposal must be made according to official regulations.

## SECTION 14. TRANSPORT INFORMATION:

**DOT:** Not regulated.

**IMDG:** Not regulated.

**IATA:** Not regulated.

## SECTION 15. REGULATORY INFORMATION:

### US Federal Regulations:

SARA 302 Components: None present or none present in regulated quantities.

SARA 313 Components: None present or none present in regulated quantities.

SARA 311/312 Hazards: No data available.

### US State Regulations:

California Prop. 65: No ingredient regulated by CA Prop 65 present.

Massachusetts RTK: No ingredient regulated by MS Right to Know Act present.

New Jersey & Community RTK: Water CAS No. 7732-18-5

Pennsylvania RTK: Water CAS No. 7732-18-5

Rhode Island RTK: No ingredient regulated by RI Right to Know Act present.

# Safety Data Sheet

Trade Name: TL Buffer

Revision Date: 4/29/2016

## SECTION 16. OTHER INFORMATION:

**Issue Date:** 9/10/2015

**Revision Date:** 4/29/2016

**Version:** 3.0

**Further Information:** No data available.

**Disclaimer:** The above information is believed to be correct but does not claim to be all inclusive. This document shall be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Omega Bio-tek, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.