

BEADS PEN

PRECISE PRESSING, STRICTLY STERILIZED

The beads Pen is an innovative product designed for use with a homogenizer. The grinding beads are pre-installed in the beads pen, and the desired number of grinding beads can be dispensed according to the requirements. It perfectly solves the problem of adding beads into the grinding tube, no need to use tweezers anymore.No need to worry about the contamination problem of enzyme-free zirconia beads after opening the cap anymore.

Product Application:

It solves the problem of adding beads into the grinding tube, dispensing the desired number of beads as needed.



>> Product Features

- Directly press to use, convenient. No need to clip, no pouring.
- Use the whole process of clean, reduce RNA, protein degradation area.
- Two materials optional: zirconia & steel beads.
- Three specifications optional: 2mm, 3mm, 4mm.
- Beads pen capacity: 20 mL.

Product Name	Cat. No.	Spec.	Beads Output Per Time / Piece	Volume	Pen Color	Beads Color
Beads Pen (Sterile, RNase & DNase free)	SYM-02Z	Piece (Zirconia Beads 2 mm)	1, 2, 3, 4, 5	20 mL 70g (2000 pcs)	Green	Dark brown, white
	SYM-03Z	Piece (Zirconia Beads 3 mm)	1, 2, 3, 4	20 mL 70g (600 pcs)	Green	Dark brown, white
	SYM-04Z	Piece (Zirconia Beads 4 mm)	1, 2, 3	20 mL 70g (400 pcs)	Green	Dark brown, white
Beads Pen (Sterile)	SYM-02B	Piece (Steel Beads 2 mm)	1, 2, 3, 4, 5	20 mL 90g (2000 pcs)	Blue	Silver
	SYM-03B	Piece (Steel Beads 3 mm)	1, 2, 3, 4	20 mL 90g (600 pcs)	Blue	Silver
	SYM-04B	Piece (Steel Beads 4 mm)	1, 2, 3	20 mL 90g (400 pcs)	Blue	Silver

The above beads pen have undergone irradiation sterilization and the zirconia beads have turned dark brown, RNase and DNase free.

GRINDING BEADS

Grinding Beads (Steel)

1 mm

2 mm

3 mm

4 mm

Grinding Beads (Glass)

0.1 mm

0.5 mm

Grinding Beads (Zirconia, RNase Free)

0.1 mm

0.5 mm

1 mm

2 mm

3 mm

4 mm

Name	Cat. No.	Spec.	Applicable Model
Grinding Beads (Steel)	G0101-200G	1 mm 200 g	Universal
	G0102-200G	2 mm 200 g	Universal
	G0103-200G	3 mm 200 g	Universal
	G0104-200G	4 mm 200 g	Universal
Grinding Beads (zirconia)	G0201-150G	1 mm 150 g	Universal
	G0202-150G	2 mm 150 g	Universal
	G0203-150G	3 mm 150 g	Universal
	G0204-150G	4 mm 150 g	Universal
	G0210-100G	0.1 mm 100 g	Universal
	G0250-100G	0.5 mm 100 g	Universal
Grinding Beads (Glass)	G0310-40G	0.1 mm 40 g	Universal
	G0350-40G	0.5 mm 40 g	Universal
Grinding Tube 2.0 mL	HT-200-M	500 pcs/bag	Universal